




## Quality Check on FFPE Samples

The CIRC team developed a QC step, using a modified RAPD method to assess FFPE quality (Figure 1). Briefly, 5 ng of DNA from FF and FFPE DNA were amplified with RAPD primers, KAPA2G Fast HotStart Taq Polymerase (Kapa Biosystems), MgCl<sub>2</sub>, and dNTPs. PCR was performed for 40 cycles.

To demonstrate concordance, the CIRC team also performed a QC check using an established, gene-specific, quantitative PCR method (Kapa Biosystems). The results of the QC step were in concordance with the QC results following the modified RAPD method.



**Figure 1: Amplification Results of QC Step**—Gel electrophoresis of RAPD-amplified DNA demonstrated that high-quality FFPE samples show a fingerprint pattern that resembles the profile of FF samples, with the presence of a prevalent, bright, 500 bp fragment as the highest amplicon size. Conversely, low-quality FFPE samples (denoted by the blue arrow) amplify only short fragments, with faint bands of less than 200 bp.

## Library Preparation and Exome Sequencing

The team used the Nextera Rapid Capture Exome (NRCE) Kit to prepare DNA libraries from FF tumor tissue, FFPE tumor tissue, and matched normal peripheral blood. FFPE samples were chosen to represent 3 high-quality samples (HQ-FFPE) and 1 low-quality sample (LQ-FFPE). Starting with 100 ng of genomic DNA, they followed the Nextera Rapid Capture protocol to tagment (tag and fragment) DNA, denature libraries to single-stranded DNA, hybridize to biotinylated probes, and pool libraries from FF and FFPE samples.

“Our experience with the NRCE Kit has been very positive,” said Dr. Astolfi. “We wanted to develop a protocol that was easy to implement and that could be performed in any biomolecular laboratory. NRCE possesses its own DNA tagmentation reaction and does not require initial DNA fragmentation by Covaris or another method. The ability to pool a different number of samples per reaction also makes NRCE flexible for users with distinct needs. It’s important to note that depending on the average size of the fragments, labs will need to determine the amount of input DNA and the concentration of the transposon enzyme to use.”

## Modified RAPD Method for FFPE QC

### PCR Mix (25 µl total volume)

2.5 µl	10X buffer
2 µl	25 mM MgCl <sub>2</sub>
0.5 µl	10 mM dNTPs
1.25 µl	10 µM forward primer (5'-AATCGGGCTG-3')
1.25 µl	10 µM reverse primer (5'-GAAACGGGTG-3')
0.1 µl	Taq polymerase

### Thermal Cycler Protocol

95 °C	2 minutes	} 40 cycles
95 °C	15 seconds	
37 °C	15 seconds	
72 °C	1 minute	
72 °C	30 seconds	

Library preparation yielded sufficient material for both FF and HQ-FFPE samples, even though average library size was smaller (Figure 2). For LQ-FFPE samples, the team pooled 2 different reactions to yield 500 ng of tagged DNA.

“Our experience with the NRCE Kit has been very positive.”

Following library preparation, the team sequenced the libraries in a 2 x 100 bp whole-exome sequencing run on the HiScanSQ™ System. Sequencing generated an average of 57 million reads for FF and 56 million reads for HQ-FFPE samples. All other sequencing parameters were comparable between FF and HQ-FFPE samples (Table 1). Conversely, LQ-FFPE failed to meet sequencing thresholds and quality parameters, with a much lower data yield (14 million reads) and low average coverage.

## Analysis and Results

### Alignment and Variant Calling

After quality trimming and adapter removal, sequence reads were mapped to the human reference genome hg19 with BWA software (Figure 3). Alignments were processed using SAMtools to remove PCR duplicates and using GATK to perform local realignment around the indel position, base quality score recalibration, and indel calling. Variant calling was performed with MuTect to identify single nucleotide variants (SNVs). Only variants that were covered at > 10x depth and mapped onto the 37 Mb NRCE target region were analyzed further. “We observed that the quality of the variant call is lower in LQ-FFPE samples,” said Valentina Indio, PhD, bioinformatics scientist at CIRC. “In contrast, the number of variant calls and the quality of base calls in HQ-FFPE samples are similar to FF samples.”



