


**Figure 3: NextSeq Series Sequencing Applications**—The flexible NextSeq Series enables researchers to transition seamlessly between sequencing applications to advance their research.

## Easy Amplicon Panel Customization With DesignStudio

The NextSeq Series Amplicon Sequencing Solution supports targeted resequencing of single exons, and fixed or custom panels. TruSight fixed panels expedite research and conserve limited budgets by assessing multiple genes simultaneously at approximately the same cost as a single gene assay. For example, hematological cancers can be studied using the TruSight Myeloid Sequencing Panel, a fixed panel that employs TSCA assay chemistry and consists of 54 genes mutated frequently in myeloid malignancies.

Researchers can also use DesignStudio to create customized panels with optimized content for their studies (Figure 4). This personalized, easy-to-use, web-based assay design tool provides dynamic feedback to enhance target region coverage, reducing the time required to design custom projects. By building a custom panel, users can focus their sequencing coverage on regions of interest and gain more confidence in detection of germ line and somatic variants. Using DesignStudio, researchers can create panels with a range of amplicon sizes optimized for different sample types (eg, formalin-fixed, paraffin-embedded) and reference genomes (eg, human, mouse, rat, and bovine).



**Figure 4: DesignStudio Enables Easy Creation of Custom Amplicon Sequencing Panels**—Building custom amplicon panels in DesignStudio enables researchers to sequence regions of interest to detect germline and somatic variants.


## Integrated Library Preparation

TSCA is used to prepare custom amplicon panels for sequencing. Following panel creation using DesignStudio, oligonucleotide probes are synthesized at Illumina and pooled into a Custom Amplicon Tube. The TSCA assay kit generates up to 1536 amplicons per reaction, incorporating integrated indexes to support up to 96 samples per run. It allows researchers to sequence hundreds of genomic regions covering as little as 2 kb or up to 650 kb of cumulative sequence.

For Research Use Only. Not for use in diagnostic procedures.

## Fast, Accurate Variant Detection

Leveraging Illumina NGS technology, the NextSeq Series delivers industry-leading sequencing accuracy of > 75% of sequenced bases over Q30<sup>†</sup> at 2 × 150 bp. NextSeq Systems can successfully sequence even the most difficult regions (GC-rich, homopolymers), yielding a higher percentage of true variants than other high-throughput desktop sequencers (Figure 5). Its low false-positive and false-negative rates drastically reduce the time and cost of downstream validation. New NextSeq v2 reagent kits are optimized to improve base calling and data quality even further.



**Figure 5: NextSeq Series Delivers High-Quality Sequencing Through Homopolymer Regions**—The NextSeq Series can sequence successfully through GC-rich and homopolymer regions to capture the full value of amplicon sequencing.

The easy workflow of the NextSeq Series simplifies amplicon sequencing. It takes less than 20 minutes to load and initiate the system. Sequencing of up to 96 samples is completed in as little as 26 hours. The flexible, scalable NextSeq Series with its dual sequencing output modes (Mid and High) and ability to handle a range of sample sizes enables researchers to tune and optimize their amplicon studies easily (Table 1).

**Table 1: NextSeq Series Supports Various Study Sizes**

Panel Size	Average Coverage Depth	Sequencing Output	Samples / Run
270 amplicons	≥ 5000x	2 × 150 (Mid Output)	96
1300 amplicons	≥ 1000x	2 × 150 (Mid Output)	96
1300 amplicons	≥ 5000x	2 × 150 (Mid Output)	24
1300 amplicons	≥ 5000x	2 × 150 (High Output)	72

Based on industry-leading Illumina sequencing by synthesis chemistry, the NextSeq Series enables researchers to compare and integrate the data it generates with data from studies performed on other Illumina systems. For example, NextSeq Series amplicon sequencing data can be integrated with data generated on the Illumina MiSeq<sup>®</sup> or HiSeq<sup>®</sup> Systems (Table 2).

<sup>†</sup> Q30 = 1 error in 1,000 base calls or an accuracy of 99.9%



Table 3: NextSeq Series Performance Parameters

Flow Cell Configuration	Read Length (bp)	Output (Gb)	Run Time	Data Quality	Required Input TruSeq Custom Amplicon Assay
High-Output Flow Cell	2 × 150	100–120	29 hours	> 75% > Q30	≥ 50 ng High-Quality Genomic DNA ≥ 250 ng FFPE Genomic DNA
Up to 400 M single reads	2 × 75	50–60	18 hours	> 80% > Q30	
Up to 800 M paired-end reads	1 × 75	25–30	11 hours	> 80% > Q30	
Mid-Output Flow Cell	2 × 150	32.5–39	26 hours	> 75% > Q30	
Up to 130 M single reads	2 × 75	16.25–19.5	15 hours	> 80% > Q30	
Up to 260 M paired-end reads	2 × 75	16.25–19.5	15 hours	> 80% > Q30	

Total times include cluster generation, sequencing, and base calling on a NextSeq System. Install specifications are based on Illumina PhiX control library at supported cluster densities (between 129 and 165 K/mm<sup>2</sup> clusters passing filter). Actual performance parameters may vary based on sample type, sample quality, and clusters passing filter. The percentage of bases > Q30 is averaged over the entire run.

Table 4: NextSeq Series Specifications

Instrument Configuration
RFID tracking for consumables
Instrument Control Computer (Internal) <sup>a</sup>
Base Unit: Dual Intel Xeon ES-2448L 1.8 GHz CPU
Memory: 96 GB RAM
Hard Drive: 750 GB
Operating System: Windows 7 embedded standard
Operating Environment
Temperature: 19°C to 25°C (22°C ± 3°C)
Humidity: Noncondensing 20%–80% relative humidity
Altitude: Less than 2,000 m (6,500 ft)
Air Quality: Pollution degree rating of II
Ventilation: Up to 2,048 BTU/hr @ 600 W
For Indoor Use Only
Light Emitting Diode (LED)
520 nm, 650 nm; Laser diode: 780 nm, Class IIIb
Dimensions
WxDxH: 53.3 cm × 63.5 cm × 58.4 cm (21.0 in × 25.0 in × 23.0 in)
Weight: 83 kg (183 lbs)
Crated Weight: 151.5 kg (334 lbs)
Power Requirements
100–120 VAC 15 A
220–240 VAC 10 A
Radio Frequency Identifier (RFID)
Frequency: 13.56 MHz
Power: Supply current 120 mA, RF output power 200 mW
Product Safety and Compliance
NRTL certified IEC 61010-1
CE marked
FCC/IC approved

a. Computer specifications are subject to change.

## Ordering Information

System Name	Catalog No.
NextSeq 500 System	SY-415-1001
NextSeq 550 System	SY-415-1002
Output Kit Name	Catalog No.
NextSeq 500 Mid-Output Kit (150 cycles)	FC-102-1001
NextSeq 500 Mid-Output Kit (300 cycles)	FC-404-1003
NextSeq 500 High-Output Kit (75 cycles)	FC-404-1005
NextSeq 500 High-Output Kit (150 cycles)	FC-404-1002
NextSeq 500 High-Output Kit (300 cycles)	FC-404-1004
NextSeq 500 Mid-Output v2 Kit (150 cycles)	FC-404-2001
NextSeq 500 High-Output v2 Kit (150 cycles)	FC-404-2002
NextSeq 500 Mid-Output v2 Kit (300 cycles)	FC-404-2003
NextSeq 500 High-Output v2 Kit (300 cycles)	FC-404-2004
NextSeq 500 High-Output v2 Kit (75 cycles)	FC-404-2005
Library Preparation Kit Name	Catalog No.
TruSeq Custom Amplicon v1.5 Kit (96 samples)	FC-130-1001
TruSeq Custom Amplicon Index Kit (96 indexes, 384 samples)	FC-130-1003
TruSeq Index Plate Fixture Kit	FC-130-1005
TruSeq Index Plate Fixture and Collar Kit (2 each)	FC-130-1007

## Reference

1. Raczy C, Petrovski R, Saunders CT, et al. Isaac: Ultra-fast whole genome secondary analysis on Illumina sequencing platforms. *Bioinformatics*. 2013;29:2041-2043.



AGAATGATAACAGTAAACACACTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCAACGTACCCTAAACGAACGATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTTCTTCTGTTAACCTTAAAGATTACTTGATCCACTGATTCAACG  
TCAACGTACCGTAAACGAACGATCATTAAAGATTACTTGATCCACTGATTCAACGTACCCTAAACGAACGATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTTCTTCTGTTAACCTTAAAGATTACTTGATCCACTGATTCAACG  
CGACGAAAAAGATGATAACAGTAAACACACTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCAACGTACCCTAAACGAACGATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTTCTTCTGTTAACCTTAAAGATTACTTGATCCACTGATTCAACG  
AAGCGTACCATTAAGAGCTACCGTCCAAACAGTAAACACTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCAACGTACCCTAAACGAACGATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTTCTTCTGTTAACCTTAAAGATTACTTGATCCACTGATTCAACG  
AGAAATGATAACAGTAAACACACTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCAACGTACCCTAAACGAACGATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTTCTTCTGTTAACCTTAAAGATTACTTGATCCACTGATTCAACG  
GATTACTTGATCCACTGATTCAACGTAAAGATTACTTGATCCACTGATTCAACGTACCCTAAACGAACGATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTTCTTCTGTTAACCTTAAAGATTACTTGATCCACTGATTCAACG  
CGTTCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTTCTTCTGTTAACCTTAAAGATTACTTGATCCACTGATTCAACGTACCCTAAACGAACGATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTTCTTCTGTTAACCTTAAAGATTACTTGATCCACTGATTCAACG