

Local Run Manager TruSight Oncology Comprehensive (JP) Analysis Module

Workflow Guide

This document and its contents are proprietary to Illumina, Inc. and its affiliates ("Illumina"), and are intended solely for the contractual use of its customer in connection with the use of the product(s) described herein and for no other purpose. This document and its contents shall not be used or distributed for any other purpose and/or otherwise communicated, disclosed, or reproduced in any way whatsoever without the prior written consent of Illumina. Illumina does not convey any license under its patent, trademark, copyright, or common-law rights nor similar rights of any third parties by this document.

The instructions in this document must be strictly and explicitly followed by qualified and properly trained personnel in order to ensure the proper and safe use of the product(s) described herein. All of the contents of this document must be fully read and understood prior to using such product(s).

FAILURE TO COMPLETELY READ AND EXPLICITLY FOLLOW ALL OF THE INSTRUCTIONS CONTAINED HEREIN MAY RESULT IN DAMAGE TO THE PRODUCT(S), INJURY TO PERSONS, INCLUDING TO USERS OR OTHERS, AND DAMAGE TO OTHER PROPERTY, AND WILL VOID ANY WARRANTY APPLICABLE TO THE PRODUCT(S).

ILLUMINA DOES NOT ASSUME ANY LIABILITY ARISING OUT OF THE IMPROPER USE OF THE PRODUCT(S) DESCRIBED HEREIN (INCLUDING PARTS THEREOF OR SOFTWARE).

© 2025 Illumina, Inc. All rights reserved.

All trademarks are the property of Illumina, Inc. or their respective owners. For specific trademark information, refer to www.illumina.com/company/legal.html.

Table of Contents

| | |
|--|-----------|
| Overview | 1 |
| About This Guide | 1 |
| Enter Run Information | 2 |
| TSO Comprehensive (JP) Analysis Module Information | 2 |
| Set Run Parameters | 3 |
| Specify Samples for the Run | 3 |
| Edit Run and Initiate Sequencing | 7 |
| Analysis Methods | 7 |
| Run Quality Control | 8 |
| FASTQ Generation | 8 |
| DNA Alignment and Error Correction | 8 |
| Small Variant Calling | 9 |
| Small Variant Annotation | 10 |
| Gene Amplification Calling | 11 |
| Tumor Mutational Burden | 11 |
| Microsatellite Instability Status | 12 |
| Quality Control for DNA Sample Libraries | 12 |
| Low Depth Reporting for DNA Sample Libraries | 12 |
| RNA Alignment | 13 |
| RNA Fusion Calling | 13 |
| RNA Splice Variant Calling | 14 |
| RNA Fusion Merging | 14 |
| RNA Splice Variant Annotation | 14 |
| Quality Control for RNA Sample Libraries | 15 |
| Transcripts | 15 |
| Control Reporting | 15 |
| Tumor Profiling of Variants | 15 |
| Analysis Output | 17 |
| Files | 18 |
| Results Reports | 18 |
| TruSight Oncology Comprehensive (JP) Report | 18 |
| Sample Sheet | 38 |
| Control Output Report | 39 |
| Metrics Output | 43 |

| | |
|---|-----------|
| Low Depth Report | 46 |
| Output Folder Structure | 47 |
| View Analysis Results | 47 |
| Samples & Results | 48 |
| Report Regeneration | 50 |
| Regenerate a Report or Requeue Analysis | 50 |
| Viewing Report Regeneration Results | 51 |
| Appendix A QC Metrics Flowchart | 52 |
| Appendix B QC Metrics | 54 |
| Quality Control Metrics | 54 |
| DNA Expanded Metrics | 58 |
| RNA Expanded Metrics | 58 |
| Appendix C TSO Comprehensive (JP) Report Reference | 60 |
| Appendix D MNVs, Indels, and Deletions in EGFR and RET Detectable by Phased Variant Caller | 61 |
| Appendix E Install a Knowledge Base | 92 |
| Appendix F Cybersecurity | 94 |
| Antivirus or Antimalware Software | 94 |
| TSO Comprehensive Assay Security Certificate | 94 |
| Regenerate Security Certificate | 95 |
| Revision History | 96 |
| Technical Assistance | 97 |

Overview

The Illumina® Local Run Manager TruSight™ Oncology Comprehensive (JP) Analysis Module analyzes sequencing reads of DNA and RNA libraries prepared using the TruSight Oncology Comprehensive (JP) (TSO Comprehensive (JP)) assay. Refer to the *TruSight Oncology Comprehensive (JP) Assay Workflow Guide (document # 200041566)* for the TSO Comprehensive (JP) assay intended use.

The TSO Comprehensive analysis module supports run setup, sequencing, analysis, and reporting for the prepared DNA and RNA libraries. For patient samples, the TSO Comprehensive analysis module generates:

- A TSO Comprehensive (JP) report for each patient sample including tumor profiling, and quality control results (available in PDF and JSON formats).
- A low depth report file in tab separated format (*.tsv) for each patient sample. The file includes a list of genomic positions (annotated with gene symbols) having insufficient sequencing depth to rule out the presence of a small variant in a DNA library.
- A quality control metrics file (*.tsv) including analysis status and quality control metrics for all patient samples in a sequencing run.

For controls, the TSO Comprehensive analysis module generates a control output report (*.tsv) including quality control results for any controls in the sequencing run.

The TSO Comprehensive analysis module comprises the TSO Comprehensive (JP) Software Suite, a Knowledge Base (KB), and a TSO Comprehensive (JP) Claims Package. The KB and the TSO Comprehensive (JP) Claims Package are installed into the TSO Comprehensive analysis module. For the analysis module part number, refer to *TruSight Oncology Comprehensive (JP) Assay Workflow Guide (document # 200041566)*.

About This Guide

This guide provides instructions for setting up run parameters for sequencing and analysis using the TSO Comprehensive analysis module. Use of the software requires basic knowledge of the current Windows operating system and web browser-based user interface. For information about the TSO Comprehensive analysis module dashboard and system settings, refer to the *NextSeq 550Dx Instrument Reference Guide for Japan (document # 1000000009513)*.

Enter Run Information

Use the TruSight Oncology Comprehensive (JP) Analysis Module software to set up TSO Comprehensive (JP) runs.

Before beginning the run, make sure that a compatible KB is installed. If a compatible KB is not installed, refer to [Appendix E Install a Knowledge Base on page 92](#).

Enter run and sample setup information directly into the TSO Comprehensive analysis module.

TSO Comprehensive (JP) Analysis Module Information

The TSO Comprehensive analysis module includes analysis module, KB, and claims package version information on the Modules & Manifests screen.

1. Open TSO Comprehensive analysis module on your instrument.
2. Use the Tools menu to navigate to the Modules & Manifests screen.
3. Select **TSO Comp (JP)**.

The Modules & Manifests screen displays the following installation information:

- **Device Identifier**—A unique device identifier for the installed TSO Comprehensive analysis module and associated Claims Package. The installed KB version does not impact this identifier.
- **Product Identifier**—The version of the installed TSO Comprehensive analysis module.
- **Modified On**—The date and time that the TSO Comprehensive analysis module itself was last installed or updated.
- **Sequencing Run Settings**—Displays the read type (paired-end) and read length settings associated with the TSO Comprehensive analysis module.
- **Claims Installed**—Displays the version of the installed claims package.
- **TSO Comprehensive Security Certificate**—HTTPS certificate specific to this instrument. Required for remote access using a web browser of this instrument from another machine in the same network. Refer to [Appendix F Cybersecurity on page 94](#) for installation instructions.
- **Knowledge Base Version**—Refer to [Appendix E Install a Knowledge Base on page 92](#) for instructions on installing or updating the KB. This section includes KB installation information for the following fields:

| Field | Description |
|----------------|---|
| Name | KB name |
| Version | KB version |
| RefSeq Version | RefSeq version included in the KB. The RefSeq version shown indicates which NCBI Homo sapiens Annotation Release ¹ it originates from. |
| Published | KB publication date |
| Installed | KB installation date |
| State | KB installation State. Displays as Ready when installation is complete. |

¹ NCBI Homo sapiens Updated Annotation Release 105.20201022.

https://www.ncbi.nlm.nih.gov/genome/annotation_euk/Homo_sapiens/105.20201022.

Set Run Parameters

1. Log in to Local Run Manager on the instrument or from a networked computer.
2. Select **Create Run**, and then select **TSO Comp (JP)**.
3. In the **Run Name** field, enter a run name that identifies the run from sequencing through analysis with the following criteria.
 - 1–40 characters.
 - Only alphanumeric characters, underscores, or dashes.
 - An alphanumeric character must precede and follow dashes or underscores.
 - Unique across all runs on the instrument.
4. **[Optional]** In the **Run Description** field, enter a run description to help identify the run with the following criteria.
 - 1–150 characters.
 - Only alphanumeric characters or spaces.
 - An alphanumeric character must precede and follow spaces.

Specify Samples for the Run

Specify samples for the run using the following options:

- **Enter samples manually**—Use the blank table at the bottom of the Create Run screen.
- **Import sample sheet**—Navigate to an external file in a comma-separated values (*.csv) format.

Enter Samples Manually

1. Enter a unique sample ID in the Sample ID field with the following criteria. **Add all controls before intended use samples.** Refer to [Controls on page 5](#) for more information.

- 1–25 characters.
 - Only alphanumeric characters, underscores, or dashes.
 - An alphanumeric character must precede and follow dashes or underscores.
2. **[Optional]** Enter a sample description in the Sample Description field with the following criteria.
 - 1–50 characters.
 - Only alphanumeric characters, dashes, underscores, or spaces.
 - An alphanumeric character must precede and follow dashes, spaces, or underscores.
 3. Select an index for the DNA library and/or RNA library prepared from the sample. Refer to the *TruSight Oncology Comprehensive (JP) Assay Workflow Guide (document # 200041566)* for number of libraries and index ID selection.
 - Make sure that RNA and DNA samples are in separate columns.
 - The DNA i7+i5 Sequence field autopopulates after selecting a DNA Index ID. The RNA i7+i5 Sequence field autopopulates after selecting an RNA Index ID.
 - For a DNA sample library, select a unique index ID (UPxx or CPxx indexes) from the DNA Index ID drop-down list.
 - For an RNA sample library, select a unique index ID (UPxx only) from the RNA Index ID drop-down list.
 - If there are three libraries in total in the run, follow the index selection guidelines in the *TruSight Oncology Comprehensive (JP) Assay Workflow Guide (document # 200041566)*.
 4. Use the Tumor Type field to assign a tumor type for each sample, selecting the most specific tumor type available. Refer to [Select a Tumor Type on page 5](#).
 5. Use the Tumor Type field to assign one of the following control types for each control. Refer to [Controls on page 5](#).
 - DNA External Control (DNA Positive Control)
 - DNA No-Template Control
 - RNA External Control (RNA Positive Control)
 - RNA No-Template Control
 6. Assign Sex. For Controls, Sex is Unknown.
 7. **[Optional]** Select **Export to CSV** to export sample information to a file.
 8. Review the information on the Create Run Screen. Incorrect information can impact results.
 9. Select **Save Run**.

Import Samples

1. Select **Import CSV** and browse to the location of the sample information file. There are two types of files that you can import.

- Select **Download CSV** on the Create Run screen to download a new sample information template. The CSV file contains the required column headings and format for import. Enter sample information in each column for the samples in the run. For the Tumor Type column, enter the tumor type term or associated code (refer to [Download Tumor Types on page 7](#)). The Tumor Type field is also used to designate samples as controls (refer to [Controls on page 5](#)).
 - Use the file of sample information that was exported from the TSO Comprehensive analysis module using the Export to CSV feature.
2. On the Create Run screen, review the imported information.
Incorrect information can impact results.
 3. **[Optional]** Select **Export to CSV** to export sample information to an external file.
 4. Select **Save Run**.

Controls

TSO Comprehensive (JP) requires the use of TruSight Oncology Controls. Designating a sample as a control automatically sets the Sex of the sample to Unknown. To designate a control sample, select one of four control types from the Tumor Type field:

- DNA External Control (positive DNA control)
- RNA External Control (positive RNA control)
- DNA No-Template Control
- RNA No-Template Control

Refer to [Select a Tumor Type on page 5](#) for more information on setting tumor types for all types of samples during run setup.

Specify one of each control type within a run. Select a DNA library for a DNA External Control or a DNA No-Template Control. Select an RNA library for an RNA External Control or an RNA No-Template Control. DNA or RNA No-Template controls are not counted against the maximum number of libraries in a run.

Refer to the *TruSight Oncology Comprehensive (JP) Assay Workflow Guide (document # 200041566)* for more information on using control samples.

Select a Tumor Type

A tumor type must be specified for each sample. Except for control types, the available tumor types are derived from the installed KB and might change with updated versions of the KB.




CAUTION

Incorrect selection of tumor type can cause incorrect results. Resolve any warnings that appear when specifying tumor types to avoid analysis failure.

The tumor type terms are part of a hierarchical disease ontology in the KB, which is constructed as a set of parent-child relationships. For example, the term non-small cell lung cancer is a child of lung cancer because non-small cell lung cancer is a type of lung cancer. [Figure 1](#) depicts a subset of an example disease ontology, showing solid tumor as the root term, and the terms associated with lung cancer and thyroid cancer (other tumor types are not shown). A term that is connected through parent-child relationships to lower-level terms is called an ancestor. The connected lower-level terms are descendants of the ancestor term. For example, lung cancer is an ancestor of adenocarcinoma of lung and small cell lung cancer, and medullary thyroid carcinoma is a descendant of both thyroid carcinoma and solid tumor.

Figure 1 Example of a Disease Ontology Subset



The selected tumor type for a patient sample impacts:

- Which tumor profiling variants are included in the TSO Comprehensive (JP) report. Refer to [Tumor Profiling of Variants on page 15](#).

Select a tumor type using the Create Run screen. The tumor type can also be set by importing a CSV file containing a tumor type (refer to [Import Samples on page 4](#)).

1. Double-click the Tumor Type cell to view the available tumor types. Available tumor types are displayed in an alphabetized hierarchical list. The Tumor Type field is also used to designate a control type for control samples (refer to [Controls on page 5](#)).
2. Use the list or search bar at the top of the Tumor Type window to select the desired tumor type.

Download Tumor Types

A full list of available tumor types in TSV format can be downloaded from the Create Run screen using the **Download Tumor Types TSV** button. The list contains the following information:

- The tumor type term visible in the user interface.
- The full path of the tumor type within the tumor type hierarchy (disease ontology).
- The code used by the TSO Comprehensive analysis module to identify the tumor type.

Edit Run and Initiate Sequencing

For instructions on editing the run information and initiating a sequencing run, refer to the *NextSeq 550Dx Instrument Reference Guide for Japan (document # 1000000009513)*. Analysis and reporting begin after a sequencing run is complete.

For storage considerations, a sequencing run can produce 40–100 GB of output. Secondary analysis of a sequencing run can produce 100–200 GB of output.

Analysis Methods

After collecting the sequencing data, the TSO Comprehensive analysis module processes it to:

- Perform quality control.
- Detect variants.
- Determine Tumor Mutational Burden (TMB) and Microsatellite Instability (MSI) status (not approved for Japan).
- Assess the clinical significance and potential clinical significance of detected variants.
- Report results.

The following sections describe the analysis methods.

Run Quality Control

Sequencing run quality metrics are evaluated to determine if they are within an acceptable range. The overall percentage of reads passing filter is compared to a minimum threshold. For Read 1 and Read 2, the average percentage of bases \geq Q30, which gives a prediction of the probability of an incorrect base call (Q-score), are also compared to a minimum threshold. If values for each of these three metrics meet the specifications, then Run QC is reported as PASS and analysis continues. If a value for any one of the metrics fails to meet the specification, then Run QC is reported as FAIL and analysis does not proceed. For more information, refer to [Quality Control Metrics on page 54](#).

FASTQ Generation

Sequencing data stored in BCL format is demultiplexed using index sequences unique to each sample added during the library preparation step to assign clusters to the library from which they originated. Each cluster contains two indexes (i5 and i7 sequences, one at each end of the library fragment). The combination of those index sequences is used to demultiplex the pooled libraries.

After demultiplexing, FASTQ files are generated. These files contain the sequencing reads for each individual sample library and the associated quality scores for each base call, excluding reads from any clusters that did not pass filter.

DNA Alignment and Error Correction

DNA alignment and error correction involve aligning sequencing reads derived from DNA sample libraries to a reference genome and correcting errors in the sequencing reads before variant calling.

The alignment step uses the Burrows-Wheeler Aligner (BWA-MEM) with the SAMtools utility to align DNA sequences in FASTQ files to the hg19 reference genome, generating BAM files (*.bam) and BAM index files (*.bam.bai).

The BAM files are further processed to remove errors (including errors introduced during PCR amplification or sequencing). Reads derived from the same unique DNA molecule are collapsed into a single representative sequence, using their unique molecular identifier (UMI) incorporated into the library fragments during library preparation.

A second round of alignment using BWA-MEM and SAMtools is performed on the UMI-collapsed reads, resulting in a second set of BAM files with corresponding BAM index files. These BAM files are used as input for gene amplification calling.

Candidate insertions and deletions are identified from the collapsed BAM alignments, and the read pairs are realigned against those candidate insertions and deletions to rescue insertions and deletions signals that may have been missed due to misalignment. Simultaneously, overlapping read pairs are stitched (bioinformatically combined) into a single consensus read. All reads are then output as a third set of BAM files with corresponding BAM index files. These BAM files are used as input for small variant calling, microsatellite instability (MSI) status determination, and DNA library quality control.

Small Variant Calling

Small variant calling is performed for DNA sample libraries (excluding DNA no-template controls) to detect small variants, including single-nucleotide variants (SNVs), multi-nucleotide variants (MNVs) up to 3 base pairs (bp) in length, and insertions and deletions up to 25 bp in length. Certain MNVs, indels (one or more nucleotides replaced by one or more nucleotides and is not an SNV or MNV), and deletions might require a phasing approach to be detected. A predefined set of MNVs, indels, and deletions are detected for the EGFR and RET genes (refer to [Appendix D MNVs, Indels, and Deletions in EGFR and RET Detectable by Phased Variant Caller on page 61](#)) using a phasing approach. The phasing approach for small variant calling is limited to only these variants. The variant calling algorithms do not differentiate between variants of somatic or germline origin.

Small Variant Detection

The error-corrected BAM files (collapsed and insertions and deletions realigned) are used as input by an initial variant calling algorithm to detect small variants. The initial variant calling step results in unfiltered genome Variant Call Format (gVCF) files. gVCF files contain reference or variant case calls for each locus targeted by the TSO Comprehensive (JP) assay.

Small Variant Filtering

Candidate variants are then filtered for recurrent (assay-specific) artifacts and artifacts from sample processing (such as deamination or oxidation). To address assay-specific artifacts, an adjusted quality score is calculated by comparing the observed variant frequency against a baseline noise distribution for the same site. This distribution was derived from profiling a set of normal samples matching the intended use population of varying qualities through the TSO Comprehensive (JP) assay. To address sample-specific artifacts, the reads supporting the variant call are stratified by error rate. Reads originating from duplex/stitched reads have the lowest error rate and reads originating from simplex (nonduplex/unstitched) reads have the highest error rate. These error rates are estimated by evaluating all loci with reported variant allele frequencies below 5%. Non-reference reads at these sites are largely due to error. True somatic events, because of their relative rarity, do not significantly impact these error rate estimates. Because these read classes, duplex/stitched and simplex, have different, sample-specific error rates, confident detection of a candidate variant may require more or fewer reads as a function of that error rate. For example, at a coverage depth of 200 reads, a variant may be confidently called with three high-quality supporting reads, or with five lower-quality supporting reads.

Candidate variants that do not have sufficient read support based on this error-aware model or that have low adjusted quality scores are tagged with a LowSupport filter flag and are considered as reference calls. If the site also has insufficient coverage for variant calling (less than 100x), the variant is tagged with a LowDP filter flag and is considered as a no-call. Variants with high prevalence in COSMIC3 have lower thresholds for each of these quality metrics compared to non-COSMIC variants. This filtering step results in filtered gVCF files.

Small Variant Phasing

A phased variant caller identifies certain MNVs, indels, and deletions in the EGFR and RET genes. The algorithm identifies variants in the EGFR and RET genes that are candidates for phasing in the filtered gVCF files from the previous step and arranges the variants into local neighborhoods. It then mines the error-corrected BAM file for any evidence that these small variants occur in the same clonal sub-populations with each other (in phase with each other). Overlapping reads are clustered in the neighborhood into a minimal set of clusters that contain the same variants. Variants are detected by examining the Concise Idiosyncratic Gapped Alignment Report (CIGAR) strings in the BAM file and comparing read sequences to the reference genome sequence.

Small Variant Merging

Finally, MNVs, indels, and deletions detected by the phased variant caller are merged into the filtered gVCF files. Only those MNVs, indels, and deletions from a predefined list of variants in the EGFR and RET genes are eligible for merging into the gVCF. Refer to [Appendix D MNVs, Indels, and Deletions in EGFR and RET Detectable by Phased Variant Caller on page 61](#). MNVs, indels, and deletions from the phased variant caller take precedence over those that may exist in the gVCF from the initial variant calling step. This step results in merged gVCF files.

Small Variant Annotation

Detected small variants are annotated using the Nirvana annotation engine with information from the RefSeq database and various population databases (COSMIC, ClinVar, dbSNP, 1000 Genomes, and gnomAD). Annotation of small variants is performed multiple times independently as described in the following sections.

Static Annotation Databases for TMB Calculation

Nirvana annotates filtered small variant calls with static (not updatable) annotation databases for use by downstream TMB calculation (refer to [Tumor Mutational Burden on page 11](#)). The gVCF from the Small Variant Phasing step is used as input (refer to [Small Variant Calling on page 9](#)). Variants detected by the phased variant caller are not used for TMB calculation.

Updatable RefSeq Database for Tumor Profiling

Nirvana annotates filtered small variant calls with an updatable RefSeq database as part of a downstream Tumor Profiling of Variants process (refer to [Tumor Profiling of Variants on page 15](#)). The updatable RefSeq database is included as part of the KB and may be updated periodically to be compatible with other KB content.

Gene Amplification Calling

Gene amplification calling is performed for DNA sample libraries (excluding DNA no-template controls). An algorithm is used to identify amplified genes and calculate the fold change value for the amplification genes targeted by TSO Comprehensive (JP). A fold change for a given gene is derived from the normalized read depth of the gene in the sample relative to the normalized read depth of diploid regions from the same sample. A fold change exceeding a gene-specific cutoff is considered a gene amplification. This analysis step results in a VCF file, summarizing gene amplification status and calculated fold change for each targeted amplification gene.

Each copy number variant is reported as a fold change on normalized read depth in a testing sample relative to the normalized read depth in diploid genomes. Given tumor purity, you can infer the ploidy of a gene in the sample from the reported fold change.

Given tumor purity X%, for a reported fold change Y, you can calculate the copy number n using the following equation:

$$n = [(200 * Y) - 2 * (100 - X)] / X$$

For example, a tumor purity at 30% and a MET with fold change of 2.2x indicates that 10 copies of MET DNA are observed.

Tumor Mutational Burden

TMB is calculated for DNA sample libraries (excluding DNA no-template controls). A TMB score is generated from the gVCF file generated by the Small Variant Filter step (refer to [Small Variant Calling on page 9](#)) and the annotations generated during Small Variant Annotations. SNVs and insertions and deletions variants are included in calculating the TMB score, which is derived from the count of non-driver somatic variants per megabase (evaluable region). Driver mutations are identified and filtered based on COSMIC count. TSO Comprehensive (JP) does not differentiate between variants of somatic or germline origin for purposes of small variant calling. Variants are flagged as likely germline for calculating the TMB score, applying a combination of population database and post-database filtering strategies. Variants that are observed frequently across population database are likely of germline origin. After database filtering, the proximity filter labels variants as germline if they are surrounded by database-labeled germline variants. Variants identified as likely germline are excluded from the TMB score calculation. The evaluable region is dynamically adjusted per sample based on sequencing depth. Genomic regions with a high background noise level are excluded from the TMB calculation. TMB is calculated as the number of somatic non-hotspot variants with VAF \geq 5% divided by the evaluable region size.

NOTE: Tertiary analysis of MSI and TMB is out of scope of approved indications in Japan. Those results are reference information.

Microsatellite Instability Status

To determine the MSI status of a sample, a total of 130 predefined MSI sites are evaluated. For each site, the repeat length distribution is compared against a panel of normal samples to see if the repeat distribution is significantly shifted. The final MSI score is calculated as the number of unstable sites divided by the total number of usable sites (sites with sufficient coverage). A sample is considered MSI-High if its MSI score is $\geq 20.00\%$ and MS-Stable if its MSI score is $< 20.00\%$.

NOTE: Tertiary analysis of MSI and TMB is out of scope of approved indications in Japan. Those results are reference information.

Quality Control for DNA Sample Libraries

DNA sample libraries (patient samples only) are assessed for potential contamination by DNA from other samples (foreign DNA) using a combination of a contamination score and a contamination p-value. In contaminated samples, there are germline variants (single nucleotide polymorphisms, or SNPs) that have VAF shifts from expected values of 0%, 50%, or 100%. The algorithm computes a log likelihood score across all common SNP positions where SNV calls are reported. The larger the contamination score, the more likely there is foreign DNA contamination. The rearrangement p-value summarizes a chromosome imbalance score, which represents the overall likelihood of the observed variant calls across each chromosome. If both the contamination score and rearrangement p-value are above predefined quality thresholds, a sample is considered to be contaminated. If contamination is detected, then DNA Library QC is reported as FAIL and no results are available for small variants, gene amplifications, MSI, and TMB. Also, a tumor profiling result is not available if it relies on DNA library QC passing.

QC metrics are used to assess the validity of small variant calling, gene amplifications, TMB, and MSI for DNA sample libraries that pass contamination quality control. If the sample library fails one or more quality metrics, then the corresponding variant type or biomarker is not reported. The associated QC category in the report header displays as FAIL. Also, a tumor profiling result may not be available if it relies on QC passing for one or more of the below QC categories.

DNA library QC results are available in the `MetricsOutput.tsv` file. Refer to [Metrics Output on page 43](#).

Low Depth Reporting for DNA Sample Libraries

A Low Depth Report is generated for each patient sample with a DNA library. The report includes a listing of genomic positions with a total sequencing depth < 100 and for which a passing small variant was not detected. These positions have insufficient sequencing depth to rule out the presence of a small variant. If there is sufficient sequencing depth of the variant allele, it is still possible to detect variants with a total sequencing depth < 100 .

Contiguous positions of low depth overlapping the same genes are combined into genomic ranges in the Low Depth Report. Each genomic range in the report is annotated with one or more RefSeq gene symbols. The RefSeq annotation is based on the RefSeq database included as part of the KB and may change with a KB update.

Refer to [Low Depth Report on page 46](#) for details on the content.

RNA Alignment

RNA alignment is performed for RNA sample libraries. RNA alignment includes preprocessing of unaligned sequencing reads, aligning sequencing reads to a reference genome, and postprocessing of aligned sequencing reads.

1. First, RNA sequences in FASTQ files are downsampled to approximately 30 million reads per RNA sample library. Downsampling is done by randomly selecting reads from the input FASTQ files following a probability distribution. Next, the ends of RNA sequences are trimmed to a maximum length of 76 base pairs.
2. Preprocessed reads are then aligned to the hg19 reference genome and candidate splice junctions are identified. This step generates BAM files and BAM index files for aligned reads, and a tab-delimited text file for candidate splice junctions.
3. Finally, duplicate reads are marked in the BAM files, such that they can be excluded from downstream steps. This step generates BAM files and BAM index files that are used as input to RNA Fusion Calling and RNA Splice Variant Calling.

RNA Fusion Calling

Fusion calling is performed for RNA sample libraries (excluding RNA no-template controls). Candidate fusions are identified from anomalous read pairs (reads aligning to different chromosomes or in unexpected orientations) in the BAM files (generated during RNA Alignment) for the fusion genes targeted by TSO Comprehensive (JP). Fusion-supporting reads are assembled into candidate fusion contigs. Candidate fusion contigs are then aligned back to the reference genome. These candidate fusion contigs are then evaluated against various filters before being reported as detected. These filters are summarized in the following table.

| Filter | Description |
|---------------|--|
| Imprecise | A low-resolution candidate, not an assembled fusion call. |
| RepeatOverlap | The fusion is tagged as overlapping with a repeat region. Only used as a filter for nonuniquely mapping fusion candidates. |
| WeakBreakend | The read/alignment evidence on one side of the fusion is weak. This filter primarily indicates that the reads only overlap the fusion by a few base pairs. Alternatively, it can indicate too much homology. |

| Filter | Description |
|------------------|--|
| DuplicateContig | The two half-contigs of the fusion are comprised of the same sequence. |
| ContigIntragenic | The realignment of half-contigs produces alignments that map to the same gene on both sides (or within 1 kb if unannotated). |
| LowQ | Unique fusion supporting reads are less than a predefined threshold (threshold is 5 for 9–16 million reads; 6 for 16–26 million reads; 7 for 26–30 million reads). |

Additional fusions may be detected through the RNA Splice Variant Calling process (refer to [RNA Splice Variant Calling on page 14](#) and [RNA Fusion Merging on page 14](#)).

RNA Splice Variant Calling

RNA splice variant calling is performed for RNA sample libraries (excluding RNA no-template controls). Candidate splice variants (junctions) from RNA Alignment are compared against a database of known transcripts and a splice variant baseline of non-tumor junctions generated from a set of normal FFPE samples from different tissue types. Any splice variants that match the database or baseline are filtered out unless they are in a set of junctions with known oncological function. If there is sufficient read support, the candidate splice variant is kept. This process also identifies candidate RNA fusions (refer to [RNA Fusion Merging on page 14](#)).

RNA Fusion Merging

Fusions identified during RNA Fusion Calling are merged with fusions from proximal genes identified during RNA Splice Variant Calling. The merged fusions are then annotated with gene symbols or names corresponding to a static database of transcripts (GENCODE Release 19). The result of this process is a set of fusion calls that are eligible for reporting.

RNA Splice Variant Annotation

Detected RNA splice variants are annotated using the Nirvana annotation engine with information from the RefSeq database. Annotation of splice variants is performed multiple times independently as described in the following sections.

Updatable RefSeq Database for Tumor Profiling

Nirvana annotates detected RNA splice variant calls with an updatable RefSeq database as part of a downstream Tumor Profiling of Variants process (refer to [Tumor Profiling of Variants on page 15](#)). Splice variants are annotated with transcript-level changes (affected exons in the gene transcript) with respect to RefSeq. The updatable RefSeq database is included as part of the KB and may be updated periodically to be compatible with other KB content.

Quality Control for RNA Sample Libraries

QC metrics are used to assess the validity of RNA sample libraries. If a QC metric is not within the acceptable range, then RNA Library QC is reported as FAIL and no results are available for fusions or splice variants. Also, a tumor profiling result is not available if it relies on RNA library QC passing.

RNA library QC results are available in the `MetricsOutput.tsv` file. Refer to [Metrics Output on page 43](#).

Transcripts

A transcript is a strand of RNA that is transcribed from DNA. That RNA can then be translated to create a protein. A gene may have multiple transcripts (for example, if different promoters are used or there are different exon splice patterns). Each transcript has a unique number. In HGVS nomenclature, a nucleotide change that affects a coding sequence can be listed with reference to a transcript. The first letter indicates the wild type allele and the second letter indicates the variant allele. For example, NM_004333.4:c.1799T>A means that at position 1799 of transcript NM_004333.4, the coding RNA encodes a T in the reference genome but is changed to an A for this variant.

Control Reporting

A control output report is generated for each analysis and includes an assessment of each control included in the run. The TSO Comprehensive analysis module automatically invalidates patient samples based on control results. Failure of DNA controls invalidates DNA samples while failure of RNA controls invalidates RNA samples.

Refer to *TruSight Oncology Comprehensive (JP) Assay Workflow Guide (document # 200041566)* for guidance on run validity and patient sample validity based on results for controls.

The control output report is available in the `ControlOutput.tsv` file. Refer to [Control Output Report on page 39](#).

Tumor Profiling of Variants

All passing, detected variants in a patient sample are matched against the installed KB to determine the genomic findings that have evidence of clinical significance or have potential clinical significance. This process is called Tumor Profiling of Variants. A genomic finding is either a single variant with evidence of clinical significance or potential clinical significance, or a grouping of variants that, when detected together, have evidence of clinical significance or potential clinical significance.

When multiple variants are listed together as a genomic finding, it means that there is evidence for clinical significance or potential clinical significance for those variants together, in at least one of the sources listed in the Informatics Details of the report. If there are multiple genomic findings, and a variant is included in more than one of these findings, then that variant may be listed more than one time on a report. A genomic finding with a single variant will only be listed at the highest level where it meets

criteria for reporting. The following examples of genomic findings and clinical meanings include multiple variants, for illustrative purposes only. For example, the genomic finding levels may vary depending on tumor type, as discussed later in this section.

- A genomic finding of a single variant NTRK1 p.(G595R) is indicated to cause resistance to one or more TRK inhibitors, in patients with a qualifying TRK fusion (FDA-approved prescribing information Larotrectinib 211710s000lbl).
- A patient in the LIBRETTO-001 clinical trial was observed to have two genomic findings, RET p.(D898_E901del) and RET p.(D903_S904delinsEP). The patient exhibited tumor response to treatment with a RET inhibitor (PMID 32846061).
- An exploratory analysis of the BOLERO-1 and -3 trials suggested that breast cancer patients with ERBB2 amplification derived clinical benefit from mTOR inhibition if the tumors displayed PI3K pathway activation or AKT1 E17K mutations (PMID 27091708).
- A BRAF p.(V600E) mutation co-occurring with TERT promoter mutation is associated with an unfavorable prognosis in papillary thyroid carcinoma per major US guidelines.

Genomic Findings with Evidence of Clinical Significance

Genomic findings with evidence of clinical significance are reported in the Genomic Findings with Evidence of Clinical Significance section of the TSO Comprehensive (JP) report (refer to [TruSight Oncology Comprehensive \(JP\) Report on page 18](#)). Genomic findings are reported in Genomic Findings with Evidence of Clinical Significance (this section) if they meet the following criteria:

- The genomic finding is associated with benefit or lack of benefit to a therapy, as evidenced by an EMA-approved drug label or FDA-approved drug label. The tumor type of the sample must be equal to or a descendant of the KB association's tumor type in the disease ontology. Refer to [Select a Tumor Type on page 5](#) for more information on the disease ontology.
- The genomic finding is associated with benefit or lack of benefit to a therapy, has diagnostic relevance, or has prognostic relevance as evidenced by published ESMO guidelines, ASCO guidelines, or other major US clinical practice guidelines. The sample's tumor type must be equal to or a descendant of the KB association's tumor type in the disease ontology. Refer to [Select a Tumor Type on page 5](#) for more information on the disease ontology.

Genomic Findings with Potential Clinical Significance

Genomic findings with potential clinical significance are reported in the Genomic Findings with Potential Clinical Significance section of the TSO Comprehensive (JP) report (refer to [TruSight Oncology Comprehensive \(JP\) Report on page 18](#)). Genomic findings are reported in Genomic Findings with Potential Clinical Significance (this section) if they meet the following criteria:

- The genomic finding meets Genomic Findings with Evidence of Clinical Significance (Level 2) criteria for genomic findings with evidence of clinical significance (for example, EMA-approved drug label, FDA-approved drug label, ESMO guidelines, ASCO guidelines, or other major US guidelines), but

only when the tumor type of the sample is not a match to the KB association's tumor type. The tumor type of the sample therefore must not be equal to and not be a descendant of the KB association's tumor type.

- The variant has a therapeutic, diagnostic, or prognostic association in the clinical literature describing a clinical study. The tumor type of the sample must be equal to or a descendant of the KB association's tumor type.
- The variant is included in eligibility criteria for an enrolling clinical trial (phase I/II, II, II/III, III, or IV) registered at clinicaltrials.gov or the EU Clinical Trials Register (EUCTR). The tumor type of the sample must be equal to or a descendant of the clinical trial's tumor type.

TMB and MSI are always reported in the Genomic Findings with Potential Clinical Significance (Level 3), regardless of the tumor type of the sample.

Leveling Changes Due to KB Updates

As clinical evidence accumulates for variants in precision oncology, KB updates are made available to reflect the changes. Variants that were initially not reportable due to lack of clinical evidence may later be reported in the Genomic Findings with Evidence of Clinical Significance or Genomic Findings with Potential Clinical Significance sections of the TSO Comprehensive (JP) report through a KB content update. Likewise, variants may move from the Genomic Findings with Evidence of Clinical Significance or Genomic Findings with Potential Clinical Significance sections of the TSO Comprehensive (JP) report or vice versa when KB content is updated. Detected variants not meeting the criteria for any level are not reported. Susceptibility or cancer risk associations are excluded from the KB and do not impact leveling. Therapeutic associations used for leveling are limited to targeted cancer therapies and immunotherapies (not including cell-based immunotherapies).

COSMIC Annotations

Variants reported in the Genomic Findings with Evidence of Clinical Significance or Genomic Findings with Potential Clinical Significance sections of the TSO Comprehensive (JP) report are annotated with a COSMIC ID, as applicable, from the Catalog of Somatic Mutations in Cancer (COSMIC) database, which is included as part of the KB.

Analysis Output

When the analysis is completed, the TSO Comprehensive analysis module generates an analysis folder in the configured output folder for the system. Refer to the *NextSeq 550Dx Instrument Reference Guide for Japan (document # 100000009513)* for more information on configuring the output folder.

To view analysis output:

1. Navigate to the directory that contains the analysis folder.
2. Open the analysis folder to view output files.

The analysis folder name is formatted as `Analysis_#` where # defaults to 1 and increments by one for each analysis requeue. A subfolder, `YYYYMMDD_HHMMSS`, is created inside the analysis folder and indicates the date and time of the analysis (for example, `20210101_145958`).

Files

This section describes the summary output files generated during analysis.

Results Reports

TSO Comprehensive (JP) reports in PDF and JSON formats are produced for each patient sample that completed analysis successfully. Results are displayed for preview on the Samples and Results tab in the Results Reports section. Samples that did not complete analysis successfully are listed with an error message. Select **Export Report** to download one TSO Comprehensive (JP) report in PDF format. Refer to the analysis output folder for TSO Comprehensive (JP) reports for all completed samples.

TruSight Oncology Comprehensive (JP) Report

The following tables describe the sections that make up the TSO Comprehensive (JP) reports produced for each patient sample in PDF and JSON formats. The PDF report is human readable, while the JSON report is built of data structures that are intended for machines to parse. Information found only in the JSON report and not reflected in the PDF report is marked as N/A for the PDF report. Variants not meeting the criteria for inclusion in the Genomic Findings with Evidence of Clinical Significance or Genomic Findings with Potential Clinical Significance (Level 2 or 3) are not included in the reports.

Refer to the *TruSight Oncology Comprehensive (JP) Assay Workflow Guide (document # 200041566)* for interpretation of results.

Refer to the JSON schema on the TSO Comprehensive (JP) support pages on the Illumina support site for additional information on the structure, fields, and possible values in the JSON report.

- **Sample, Run, and Analysis Information**—Contains general information about the patient sample and the report.

Table 1 Sample, Run, and Analysis Information

| Field in PDF report | Field in JSON report | Description |
|---------------------|------------------------------|---|
| Report Date | reportDate | Date that the report was generated. |
| N/A | reportTime | Time that the report was generated. |
| Sample ID | sampleInformation / sampleId | Sample Identifier. Patient demographics are not included. |

| Field in PDF report | Field in JSON report | Description |
|---------------------|---------------------------------------|---|
| Tumor Type | sampleInformation / tumorType | Tumor type associated with the patient sample. |
| N/A | sampleInformation / tumorTypeCode | Tumor type code associated with the patient sample. |
| N/A | sampleInformation / tumorTypePath | Tumor type path (with respect to the disease ontology) associated with the patient sample. |
| N/A | sampleInformation / tumorTypeCodePath | Tumor type code path (with respect to the disease ontology) associated with the patient sample. |
| Sex | sampleInformation / sex | Patient sex (Male, Female, or Unknown). |
| Analysis Date | sampleInformation / analysisDate | Date that the secondary analysis was completed. |
| N/A | sampleInformation / analysisTime | Time that the secondary analysis was completed. |
| Run ID | sampleInformation / analysisRunId | Sequencing run ID. |
| N/A | sampleInformation / analysisRunName | Sequencing run name. |

- **Quality Control**—Contains quality control information. For more information on how quality control is evaluated, refer to [Appendix A QC Metrics Flowchart on page 52](#).

Table 2 Quality Control

| Field in PDF report | Field in JSON report | Description |
|----------------------------|---|---|
| Run QC | qualityControl / status / (array item having label = "Run QC") | <p>Run QC (PASS, FAIL, or N/A) applies to all samples contained in a single sequencing run.</p> <ul style="list-style-type: none"> • PASS—The run is valid. • FAIL or N/A—The run is invalid. All RNA and DNA sample-specific QC statuses are N/A (DNA Library QC, DNA MSI QC, DNA Small Variant & TMB QC, DNA Copy Number Variant QC, DNA External Control & NTC, RNA External Control & NTC, and RNA Library QC) and there are no variants or biomarkers listed in the report. Refer to the <i>TruSight Oncology Comprehensive (JP) Assay Workflow Guide (document # 200041566)</i> for guidance on run validity and patient sample validity based on results for controls. |
| RNA External Control & NTC | qualityControl/status/ (array item having label = "RNA External Control & NTC") | <p>RNA External Control & NTC Control Results (PASS, FAIL, or N/A) applies to the RNA library that was sequenced.</p> <ul style="list-style-type: none"> • PASS—both the RNA External Control and RNA No-Template Control have a result of PASS. • FAIL—the RNA External Control and/or RNA No-Template Control has a result of FAIL. • N/A—the RNA library for the sample was not sequenced or reanalyzed during requeue, the RNA External Control or NTC for the run was not sequenced or reanalyzed during requeue, the sample's RNA library, external control and NTC were excluded on LRM during requeue, the run folder was manually deleted before the analysis started (in this case, all QC metrics are N/A), or Run QC had a value of FAIL. <p>If the value is FAIL or N/A, there are no RNA variant types (fusion or splice variants) in the report.</p> |

| Field in PDF report | Field in JSON report | Description |
|----------------------------|---|--|
| RNA Library QC | qualityControl / status / (array item having label = "RNA Library QC") | <p>RNA Library QC (PASS, FAIL, or N/A) applies to the RNA library that was sequenced.</p> <ul style="list-style-type: none"> PASS—the RNA library passed all the RNA-specific QC metrics. FAIL—the RNA library failed one or more of the RNA-specific QC metrics. The RNA sample read is considered zero due to an index sequence selection error in the RNA sample. N/A—the RNA library for the sample was not sequenced, the RNA library for the sample was excluded on LRM during requeue, the run folder was manually deleted before the analysis started (in this case, all QC metrics are N/A), or Run QC had a value of FAIL. <p>If the value is FAIL or N/A, there are no RNA variant types (fusion or splice variants) in the report.</p> |
| DNA External Control & NTC | qualityControl/status/ (array item having label = "DNA External Control & NTC") | <p>DNA External Control & NTC Control Results (PASS, FAIL, or N/A) applies to the DNA library that was sequenced.</p> <ul style="list-style-type: none"> PASS—both the DNA External Control and DNA No-Template Control have a result of PASS. FAIL—the DNA External Control and/or DNA No-Template Control has a result of FAIL. N/A—the DNA library for the sample was not sequenced or reanalyzed during requeue, the DNA External Control or NTC for the run was not sequenced or reanalyzed during requeue, the run folder was manually deleted before the analysis started (in this case, all QC metrics are N/A), or Run QC had a value of FAIL. <p>If the value is FAIL or N/A, there are no DNA variant types (small variants, gene amplifications) or DNA biomarkers (TMB, MSI) in the report.</p> |

| Field in PDF report | Field in JSON report | Description |
|------------------------------|--|--|
| DNA Library QC | qualityControl / status / (array item having label = "DNA Library QC") | <p>DNA Library QC (PASS, FAIL, or N/A) applies to the DNA library that was sequenced.</p> <ul style="list-style-type: none"> • PASS—the DNA library passed the contamination QC metric. • FAIL—the DNA library failed the contamination QC metric. • N/A—the DNA library for the sample was not sequenced, the run folder was manually deleted before the analysis started (in this case, all QC metrics are N/A), or Run QC had a value of FAIL. <p>If the value is FAIL or N/A, no DNA variant types (small variants, copy number variants) or DNA biomarkers (TMB, MSI) are reported.</p> |
| DNA MSI QC | qualityControl / status / (array item having label = "DNA MSI QC") | <p>DNA MSI QC (PASS, FAIL, or N/A) applies to the Solid-FFPE DNA library that was sequenced.</p> <ul style="list-style-type: none"> • PASS—the DNA library passed the MSI-specific QC metric and upstream DNA Library QC metric. • FAIL—the DNA library failed the MSI-specific QC metric. • N/A—the DNA library for the sample was not sequenced, DNA Library QC for the sample was FAIL, or Run QC had a value of FAIL. <p>If the value is FAIL or N/A, the biomarker MSI is not reported and listed as Not evaluable.</p> |
| DNA Small Variant and TMB QC | qualityControl / status / (array item having label = "DNA Small Variant & TMB QC") | <p>DNA Small Variant and TMB QC (PASS, FAIL, or N/A) apply to the DNA library that was sequenced.</p> <ul style="list-style-type: none"> • PASS—the DNA library passed the Small Variant and TMB specific QC metrics and upstream DNA Library QC metric. • FAIL—the DNA library failed one or more of the Small Variant and TMB-specific QC metrics. • N/A—the DNA library for the sample was not sequenced, DNA Library QC for the sample was FAIL, or Run QC had a value of FAIL. <p>If the value is FAIL or N/A, there are no small variants in the report, and the biomarker TMB is listed as Not evaluable.</p> |

| Field in PDF report | Field in JSON report | Description |
|----------------------------|--|--|
| DNA Copy Number Variant QC | qualityControl / status / (array item having label = "DNA Copy Number Variant QC") | <p>DNA Copy Number Variant (CNV) QC (PASS, FAIL, or N/A) applies to the DNA Solid-FFPE library that was sequenced.</p> <ul style="list-style-type: none"> • PASS—the DNA library passed all the Copy Number Variant specific QC metrics and upstream DNA Library QC metric. • FAIL—the DNA library failed one or more of the Copy Number Variant specific QC metrics. • N/A—the DNA library for the sample was not sequenced, DNA Library QC for the sample was FAIL, or Run QC had a value of FAIL. <p>If the value is FAIL or N/A, there are no gene amplifications in the report.</p> |

- **TruSight Oncology Comprehensive (JP) Analysis Module and Knowledge Base Configuration**—Contains information on the software and KB versions used when the report was generated.

Table 3 Analysis Module and KB Configuration

| Field in PDF report | Field in JSON report | Description |
|-------------------------------|--|---|
| Knowledge Base Version | softwareConfiguration / knowledgeBaseVersion | Version of the Knowledge Base installed with the TSO Comprehensive analysis module. |
| Knowledge Base Published Date | softwareConfiguration / knowledgeBasePublishedDate | Date associated with the Knowledge Base that was used to generate the report. |
| Module Version | softwareConfiguration / moduleSoftwareVersion | Version of the TSO Comprehensive analysis module used to generate the report. |
| Claims Package Version | softwareConfiguration / claimsPackageVersion | Version of the Claims Package installed with the TSO Comprehensive analysis module. |

- **Other Alterations and Biomarkers Identified**—The following two sections contain tumor profiling information for detected variants categorized into Genomic Findings with Evidence of Clinical Significance, TMB, MSI, and detected variants categorized into Genomic Findings with Potential Clinical Significance. Refer to [Tumor Profiling of Variants on page 15](#) for details on how a level is determined for detected variants.

- **Genomic Findings with Evidence of Clinical Significance**—Each entry in this section is a genomic finding, which is either a single variant with evidence of clinical significance or a grouping of variants that when detected together have evidence of clinical significance. If no variants are detected, the report displays a No Detected Variants message.

Table 4 Genomic Findings with Evidence of Clinical Significance

| Field in PDF report | Field in JSON report | Description |
|---------------------|--|--|
| Detected Variants | reportFindings / otherFindings / genomicFindingsWithEvidenceOfClinicalSignificance / results / genomicFindings / (array item for genomic finding) / variants | <p>A list of detected variants that are part of the genomic finding.</p> <p>For small variants, includes the gene symbol and protein change, transcript change, or genomic change in Human Genome Variation Society (HGVS) format, for example, NRAS p.(Q61R).</p> <p>For gene amplifications, includes the gene symbol followed by Gain, for example, ERBB2 Gain.</p> <p>For fusions, includes the symbols or names of both partner genes (from GENCODE Release 19), separated by a - or /. When separated by a -, the reported gene order corresponds to the transcribed orientation (5' to 3'). When separated by a /, orientation could not be determined. If multiple genes are overlapping a breakpoint, all are listed and delimited by semicolons.</p> <p>For splice variants, includes the gene symbol and affected exons (as applicable), for example, EGFR Exon(s) 2-7 skipped.</p> |

| Field in PDF report | Field in JSON report | Description |
|---------------------|--|--|
| Details | reportFindings / otherFindings / genomicFindingsWithEvidenceOfClinicalSignificance / results / genomicFindings / (array item for genomic finding) / variants / (array item for variant in genomic finding) | Contains a list of variant details. In the PDF report, the order of variant details corresponds to the order of variants listed for Detected Variants/Biomarkers field. Refer to Small Variant Details in Report on page 30 , Gene Amplification Details in Report on page 33 , Fusion Details in Report on page 34 , and Splice Variant Details in Report on page 36 for a list of variant detail fields. |

- Genomic Findings with Potential Clinical Significance**—TMB and MSI are both reported in this section when there is a sequenced DNA library for the sample. TMB and MSI are included as reference information and not for the purpose of diagnosis use. Each other entry in this section is a genomic finding, which is either a single variant with potential clinical significance or a grouping of variants that when detected together have potential clinical significance. If no variants are detected, the report displays a No Detected Variants message.

Table 5 Genomic Findings with Potential Clinical Significance

| Field in PDF report | Field in JSON report | Description |
|---------------------|---|--|
| TMB | reportFindings / otherFindings / biomarkers / tumorMutationalBurden | <p>TMB is a measurement of the number of estimated somatic mutations carried by tumor cells per megabase in the coding region. TMB is reported as Not evaluable if it could not be evaluated due to a QC failure or a DNA library for the sample was not sequenced.</p> <p>TMB is always included in the Genomic Findings with Potential Clinical Significance section.</p> |
| MSI | reportFindings / otherFindings / biomarkers / microsatelliteInstability | <p>MSI status. Possible values include the following:</p> <p>MS-Stable—Microsatellite stable.</p> <p>MSI-High—Microsatellite instability-high.</p> <p>Not evaluable—MSI status could not be evaluated due to a QC failure or a DNA library for the sample was not sequenced.</p> <p>MSI is always included in Genomic Findings with Potential Clinical Significance ().</p> |

| Field in PDF report | Field in JSON report | Description |
|---------------------|--|--|
| Detected Variants | reportFindings / otherFindings / genomicFindingsWithPotentialClinicalSignificance / results / genomicFindings / (array item for genomic finding) / variants / (all array items) / detectedVariantLabel | <p>A list of detected variants that are part of the genomic finding.</p> <p>For small variants, includes the gene symbol and protein change, transcript change, or genomic change in Human Genome Variation Society (HGVS) format, for example, NRAS p.(Q61R).</p> <p>For gene amplifications, includes the gene symbol followed by Gain, for example, ERBB2 Gain.</p> <p>For fusions, includes the symbols or names of both partner genes (from GENCODE Release 19), separated by a - or /. When separated by a -, the reported gene order corresponds to the transcribed orientation (5' to 3'). When separated by a /, orientation could not be determined. If multiple genes are overlapping a breakpoint, all are listed and delimited by semicolons.</p> <p>For splice variants, includes the gene symbol and affected exons (as applicable), for example, EGFR Exon(s) 2-7 skipped.</p> |

| Field in PDF report | Field in JSON report | Description |
|---------------------|---|--|
| Details | reportFindings / otherFindings / genomicFindingsWithPotentialClinicalSignificance / results / genomicFindings / (array item for genomic finding) / variants | Contains a list of variant details. In the PDF report, the order of variant details corresponds to the order of variants listed for Detected Variants/Biomarkers field. Refer to Small Variant Details in Report on page 30 , Gene Amplification Details in Report on page 33 , Fusion Details in Report on page 34 , and Splice Variant Details in Report on page 36 for a list of variant detail fields. |

- **About the Test, Informatics Details, Limitations**—Contains general information about the test and a list of limitations.

Table 6 About the Test, Informatics Details, Limitations

| Field in PDF report | Field in JSON report | Description |
|---------------------|--|---|
| About the Test | about / description | Test description. |
| Informatics Details | details / (one JSON property per subsection) | A brief description of the report sections and other informatics details. |
| Limitations | limitations / description | List of assay and report limitations. |

- **TruSight Oncology Comprehensive (JP) Gene Panel**—Contains information about the gene panel.

Table 7 TruSight Oncology Comprehensive (JP) Gene Panel

| Field in PDF report | Field in JSON report | Description |
|---------------------|--|--|
| Gene Panel | genePanel / geneList / genes / genePanel / geneList / genes / variants | The list of genes that are part of the panel, including a footnote indicating which variant types are evaluated for which genes. Small variants are called in all genes. |

- **Details in Report**—Contains information about small variants, gene amplifications, fusion variants, and splice variants.

Table 8 Small Variant Details in Report

| Field in PDF report | Field in JSON report (relative path in variant JSON object) | Description |
|---------------------|---|---|
| Type | type / value | The detailed type of variant. Possible values for small variants include: SNV —Single nucleotide variant. Insertion —Addition of nucleotides of up to 25 bp. Deletion —Removal of nucleotides of up to 25 bp. MNV —Multi-nucleotide variant, being a substitution of two or three nucleotides with the same number of nucleotides. Indel —One or more nucleotides replaced by one or more nucleotides and is not an SNV or MNV. This is commonly referred to as delins. |
| VAF | additionalInfo / (array item having label property = "VAF") / value | Variant allele frequency (as a percentage). |
| Consequence | additionalInfo / (array item having label property = "Consequence") / value | Variant consequence from the Sequence Ontology. |
| Protein Change | additionalInfo / (array item having label property = "Protein Change") / value | Change to the protein reference sequence in HGVS nomenclature, as applicable. |
| Nucleotide Change | additionalInfo / (array item having label property = "Nucleotide Change") / value | Change to the coding DNA reference sequence (RefSeq transcript) in HGVS nomenclature. If the variant does not impact a transcript, the change to the genomic reference sequence in HGVS nomenclature is included. |
| Genomic Position | additionalInfo / (array item having label property = "Genomic Position") / value | Genomic position (hg19) in chromosome:position format. Refers to the position of the first base in the reference allele. |
| Reference Allele | additionalInfo / (array item having label property = "Reference Allele") / value | Reference allele. |

| Field in PDF report | Field in JSON report (relative path in variant JSON object) | Description |
|---------------------|---|---|
| Alternate Allele | additionalInfo / (array item having label property = "Alternate Allele") / value | Alternate allele. |
| N/A | cosmicIds | List of genomic mutation IDs associated with the variant from the Catalogue of Somatic Mutations In Cancer (COSMIC) database, as applicable. |
| N/A | detailedSmallVariantData / vcfChromosome | Chromosome. |
| N/A | detailedSmallVariantData / vcfPosition | Genomic position (hg19). Refers to the position of the first base in the reference allele (detailedSmallVariantData / referenceAllele field). |
| N/A | detailedSmallVariantData / vcfRefAllele | The reference allele. |
| N/A | detailedSmallVariantData / vcfVariantFrequency | Variant allele frequency. |
| N/A | detailedSmallVariantData / annotation / transcripts | Detailed transcript-level annotations for a transcript (as applicable). Only a single preferred transcript is included. |
| N/A | detailedSmallVariantData / annotation / transcripts / (first array item) / transcript | Transcript ID. |
| N/A | detailedSmallVariantData / annotation / transcripts / (first array item) / source | Transcript source (for example, RefSeq). |
| N/A | detailedSmallVariantData / annotation / transcripts / (first array item) / bioType | An Ensembl biotype classification for the transcript. |
| N/A | detailedSmallVariantData / annotation / transcripts / (first array item) / aminoAcids | The change in amino acids, as applicable (for example, G/D). |

| Field in PDF report | Field in JSON report (relative path in variant JSON object) | Description |
|---------------------|--|--|
| N/A | detailedSmallVariantData / annotation / transcripts / (first array item) / cdnaPos | cDNA position. |
| N/A | detailedSmallVariantData / annotation / transcripts / (first array item) / codons | Codon sequence change (for example, gGt/gAt), as applicable. |
| N/A | detailedSmallVariantData / annotation / transcripts / (first array item) / cdsPos | Coding sequence position, as applicable. |
| N/A | detailedSmallVariantData / annotation / transcripts / (first array item) / exons | The exons affected by the variant, and total exon count, as applicable. For example, 4-6/7 would indicate that exons 4, 5, and 6 were affected and that this transcript contains 7 exons in total. |
| N/A | detailedSmallVariantData / annotation / transcripts / (first array item) / introns | The introns affected by the variant, as applicable. |
| N/A | detailedSmallVariantData / annotation / transcripts / (first array item) / geneld | National Center for Biotechnology Information (NCBI) gene ID. |
| N/A | detailedSmallVariantData / annotation / transcripts / (first array item) / hgnc | HUGO Gene Nomenclature Committee (HGNC) gene symbol. |
| N/A | detailedSmallVariantData / annotation / transcripts / (first array item) / consequence | Array of variant consequences from the Sequence Ontology. |
| N/A | detailedSmallVariantData / annotation / transcripts / (first array item) / hgvsC | Change to the coding DNA reference sequence (RefSeq transcript) in HGVS nomenclature, as applicable. |
| N/A | detailedSmallVariantData / annotation / transcripts / (first array item) / hgvsP | Change to the protein sequence in HGVS nomenclature, as applicable. |

| Field in PDF report | Field in JSON report (relative path in variant JSON object) | Description |
|---------------------|--|---|
| N/A | detailedSmallVariantData / annotation / transcripts / (first array item) / isCanonical | Displays true if this transcript is considered the canonical transcript of the gene, otherwise false. A canonical transcript for a gene is determined as follows: Only NM & NR transcripts are included. Transcripts for a gene are sorted in the following order: <ul style="list-style-type: none"> • Locus Reference Genomic (LRG) entries come before non-LRG entries. • Descending CDS length. • Descending transcript length. • Accession number. With this sorting, the first transcript is considered canonical. |
| N/A | detailedSmallVariantData / annotation / transcripts / (first array item) / proteinId | Protein ID. |
| N/A | detailedSmallVariantData / annotation / transcripts / (first array item) / proteinPos | Protein position. |

Table 9 Gene Amplification Details in Report

| Field in PDF report | Field in JSON report (relative path in variant JSON object) | Description |
|---------------------|---|---|
| Type | type / value | The detailed type of variant. Possible values for gene amplifications include: CNV —Copy number variant (gene amplifications are the only copy number variants listed in the report). |
| Fold Change | detailedCopyNumberVariantData / foldChange | The fold-change of normalized read depth in the sample relative to the normalized read depth in diploid genomes. |

| Field in PDF report | Field in JSON report (relative path in variant JSON object) | Description |
|---------------------|---|---|
| N/A | detailedCopyNumberVariantData / copyNumberType | Value is <DUP> for all gene amplifications. |
| N/A | detailedCopyNumberVariantData / gene | Gene symbol. |
| N/A | detailedCopyNumberVariantData / chromosome | Chromosome of the gene. |
| N/A | detailedCopyNumberVariantData / startPosition | Start position (hg19) of the gene. |
| N/A | detailedCopyNumberVariantData / endPosition | End position (hg19) of the gene. |

Annotations (positional information, consequences, etc.) provided in [Table 10](#) are based on variants that have been left-aligned to the genome in accordance with next-generation sequencing norms. The one exception to this rule is that HGVS notation is right-aligned with the respective reference sequence according to the HGVS standard. When insertions and deletions occur in low complexity genomic regions, the left-aligned and right-aligned representations might refer to different locations.

Table 10 Fusion Details in Report

| Field in PDF report | Field in JSON report (relative path in variant JSON object) | Description |
|-------------------------|---|---|
| Type | type / value | The detailed type of variant. Possible values for fusions include: Fusion |
| Breakpoint 1 | additionalInfo / (array item having label property = "Breakpoint 1") | Observed fusion breakpoint 1 in RNA. Chromosome:position format (hg19). |
| Breakpoint 2 | additionalInfo / (array item having label property = "Breakpoint 2") | Observed fusion breakpoint 2 in RNA. Chromosome:position format (hg19). |
| Fusion Supporting Reads | additionalInfo / (array item having label property = "Fusion Supporting Reads") | Count of fusion supporting reads. |

| Field in PDF report | Field in JSON report (relative path in variant JSON object) | Description |
|---------------------|---|--|
| N/A | detailedGeneFusionData / fusionDirectionalityKnownAndIndicatedByGeneOrder | Displays true when gene/breakpoint order corresponds to the transcribed orientation (5' to 3'). Displays false when orientation could not be determined. |
| N/A | detailedGeneFusionData / fusionSupportingReads | Count of fusion supporting reads. |
| N/A | detailedGeneFusionData / partner1 / gene | Symbols or name (from GENCODE Release 19) of genes overlapping Breakpoint 1. Multiple genes overlapping the same breakpoint are delimited by semicolons. |
| N/A | detailedGeneFusionData / partner1 / chromosome | Chromosome of breakpoint 1. |
| N/A | detailedGeneFusionData / partner1 / position | Position (hg19) of breakpoint 1. |
| N/A | detailedGeneFusionData / partner2 / gene | Symbols or name (from GENCODE Release 19) of genes overlapping Breakpoint 2. Multiple genes overlapping the same breakpoint are delimited by semicolons. |
| N/A | detailedGeneFusionData / partner2 / chromosome | Chromosome of breakpoint 2. |
| N/A | detailedGeneFusionData / partner2 / position | Position (hg19) of breakpoint 2. |

Table 11 Splice Variant Details in Report

| Field in PDF report | Field in JSON report (relative path in variant JSON object) | Description |
|-------------------------|---|--|
| Type | type / value | The detailed type of variant. Possible values for splice variants include: Splice Variant |
| Affected Exons | additionalInfo / (array item having label property = "Affected Exons") | The exons affected by the splice variant, as applicable. For example, 4–6 would indicate that exons 4, 5, and 6 were affected. |
| Affected Introns | additionalInfo / (array item having label property = "Affected Introns") | The introns affected by the splice variant, as applicable. For example, 3 would indicate that intron 3 was affected. |
| Transcript | additionalInfo / (array item having label property = "Transcript") | Transcript ID (RefSeq). |
| Breakpoint Start | additionalInfo / (array item having label property = "Breakpoint Start") | Observed splice variant breakpoint start in RNA. Chromosome:position format (hg19). |
| Breakpoint End | additionalInfo / (array item having label property = "Breakpoint End") | Observed splice variant breakpoint end in RNA. Chromosome:position format (hg19). |
| Splice Supporting Reads | additionalInfo / (array item having label property = "Splice Supporting Reads") | Count of splice supporting reads. |
| N/A | detailedSpliceVariantData / breakPointStartChromosome | Chromosome of breakpoint start. |
| N/A | detailedSpliceVariantData / breakPointStartPosition | Position (hg19) of breakpoint start. |
| N/A | detailedSpliceVariantData / breakPointEndChromosome | Chromosome of breakpoint end. |
| N/A | detailedSpliceVariantData / breakPointEndPosition | Position (hg19) of breakpoint end. |
| N/A | detailedSpliceVariantData / spliceSupportingReads | Count of splice supporting reads. |
| N/A | detailedSpliceVariantData / annotation / source | Transcript source (for example, RefSeq). |

| Field in PDF report | Field in JSON report (relative path in variant JSON object) | Description |
|---------------------|---|---|
| N/A | detailedSpliceVariantData / annotation / gene | Gene symbol. |
| N/A | detailedSpliceVariantData / annotation / affectedExons | The exons affected by the splice variant, and total exon count, as applicable. For example, 4-6/7 would indicate that exons 4, 5, and 6 were affected and that this transcript contains 7 exons in total. |
| N/A | detailedSpliceVariantData / annotation / affectedIntrons | The introns affected by the splice variant, and total intron count, as applicable. For example, 3/6 would indicate that intron 3 was affected and that this transcript contains 6 introns in total. |
| N/A | detailedSpliceVariantData / annotation / transcript | Transcript ID. |

Sample Sheet

File name: `SampleSheet.csv`

For each analysis, the TSO Comprehensive analysis module creates a comma-delimited sample sheet (`SampleSheet.csv`). This file contains sample information provided to the software during the run setup. These sample sheets contain a header with information about the run and descriptors for the sample libraries processed in a particular flow cell (one data row per sample library).



CAUTION

Modifying the sample sheet file causes adverse effects downstream, including incorrect results or analysis failure.

Table 12 Sample Sheet Details

| Column Name | Description |
|--------------------|--|
| Sample_ID | Sample ID with <code>-DNA</code> appended for DNA libraries or <code>-RNA</code> appended for RNA libraries. |
| I7_Index_ID | i7 index name. Refer to <i>Illumina Adapter Sequences (document # 1000000002694)</i> for details on how the sample sheet index ID maps to the index ID entered during run setup. |
| index | i7 index sequence. |
| I5_Index_ID | i5 index name. Refer to <i>Illumina Adapter Sequences (document # 1000000002694)</i> for details on how the sample sheet index ID maps to the index ID entered during run setup. |
| index2 | i5 index sequence. |
| Sample_Type | DNA or RNA. |
| Pair_ID | Sample ID (same ID is used for a DNA library and RNA library from the same sample). |
| Sample_Description | Sample description. |
| Tumor_Type | Tumor type for patient samples. Control type for controls. |
| Sex | Sex (Male, Female, or Unknown). |

Control Output Report

File name: `ControlOutput.tsv`

The control output report is a tab-delimited file that provides quality control information for any controls that were included in the run.

- If controls in the run fail, the TruSight Oncology Report includes a message that one or more controls failed and to review the control output file for details.
- If any DNA External Control or DNA No-Template Control fails, no DNA variants (small variants, gene amplifications) or biomarkers (TMB, MSI) are reported.
- If any RNA External Control or RNA No-Template Control fails, no RNA variants (fusions, splice variants) are reported.

Refer to *TruSight Oncology Comprehensive (JP) Assay Workflow Guide (document # 200041566)* for guidance on run validity and patient sample validity based on results for controls.

The control output report contains the following sections and their associated fields (run ID is included before the first section):

- **Control Types**—Contains information about each control included in the run.

Table 13 Control Types

| Field | Description |
|-----------------------|--|
| Control Type | The control type of the control. Possible values include: <ul style="list-style-type: none"> • DNA External Control • DNA No-Template Control • RNA External Control • RNA No-Template Control. |
| Sample_ID | Sample ID of the control. Value is (Not Run) if this control type was not included in the run. |
| AnalysisComplete | Indication of whether analysis completed for this control. Possible values include TRUE, FALSE, not applicable. |
| Overall Result | The QC result for the control. Possible values include PASS, FAIL, N/A. |
| Sensitivity Value | The calculated sensitivity value for the control. Represents the ratio of detected control variants to the total number of expected control variants in the control. Only applicable for the following control types: <ul style="list-style-type: none"> • DNA External Control • RNA External Control |
| Sensitivity Threshold | The minimum sensitivity value required for the control to have a QC result of PASS. Only applicable for the following control types: <ul style="list-style-type: none"> • DNA External Control • RNA External Control |

- **Analysis Details**—Contains information on the analysis.

Table 14 Analysis Details

| Field | Description |
|------------------------|---|
| Report Date | The date the control report was generated. |
| Report Time | The time the control report was generated. |
| Module Version | The version of the TSO Comprehensive analysis module. |
| Pipeline Version | The version of the analysis pipeline/workflow. |
| Claims Package Version | The version of the claims package installed with the TSO Comprehensive analysis module. |

- **Sequencing Run Details**—Contains information on the sequencing run.

Table 15 Sequencing Run Details

| Field | Description |
|-------------------------------------|--|
| Run Name | The name of the sequencing run. |
| Run Date | The date of the sequencing run. |
| Instrument ID | The unique ID associated with the sequencing instrument. |
| Instrument Control Software Version | NextSeq Control Software (NCS) version in use for the run. |
| Instrument Type | The sequencing instrument type. |
| RTA Version | Real-Time Analysis (RTA) software version in use for the sequencing run. |
| Reagent Cartridge Lot Number | The lot number of the reagent cartridge used for the run. |

- **Analysis Status**—Contains information on whether analysis completed for each control and whether any samples failed due to a software error.

Table 16 Analysis Status

| Field | Description |
|---------------------|--|
| Sample_ID | Sample ID of the control. Value is (Not Run) for control types not included in the run. |
| COMPLETED_ALL_STEPS | Indicates whether the control completed all steps of the analysis. Possible values include TRUE, FALSE, N/A. If the value is FALSE, contact Illumina technical support for more information. |

| Field | Description |
|--------------------|--|
| FAILED_STEPS | A list of any failed analysis steps due to a software error. Contact Illumina technical support for more information if any step is listed here. |
| STEPS_NOT_EXECUTED | A list of any analysis steps not executed due to a software error. Contact Illumina technical support for more information if any step is listed here. |

- **Small Variants Truth Table Results**—Contains information on the control DNA small variants in the DNA External Control (positive DNA control) that were detected or not detected (one row per control variant). N/A values are listed if the DNA External Control was not included in the sequencing run.

Table 17 Small Variants Truth Table Results

| Field | Description |
|--------------------|--|
| Detected | Indicates whether the control DNA small variant was detected in the control. Possible values include TRUE, FALSE, N/A. |
| HGNC Gene Name | HUGO Gene Nomenclature Committee (HGNC) gene symbol associated with the control DNA small variant. |
| Chromosome | Chromosome of the control DNA small variant. |
| Position | Position (hg19) of the control DNA small variant. |
| Reference Allele | Reference allele of the control DNA small variant. |
| Alternative Allele | Alternate/alternative allele of the control DNA small variant. |

- **Splice Variants Truth Table Results**—Contains information on the control RNA splice variants in the RNA External Control that were detected or not detected (one row per control variant). N/A values are listed if the RNA External Control was not included in the sequencing run.

Table 18 Splice Variants Truth Table Results

| Field | Description |
|----------------|---|
| Detected | Indicates whether the control RNA splice variant was detected in the control. Possible values include TRUE, FALSE, N/A. |
| HGNC Gene Name | HGNC gene symbol associated with the control RNA splice variant. |
| Breakpoint 1 | Chromosome and position (hg19) of the first breakpoint of the control RNA splice variant. |
| Breakpoint 2 | Chromosome and position (hg19) of the second breakpoint of the control RNA splice variant. |

- **Fusions Truth Table Results**—Contains information on the control RNA fusion variants in the RNA External Control that were detected or not detected (one row per control variant). N/A values are listed if the RNA External Control was not included in the sequencing run.

Table 19 Fusions Truth Table Results

| Field | Description |
|------------------|---|
| Detected | Indicates whether the control RNA fusion variant was detected in the control. Possible values include TRUE, FALSE, N/A. |
| HGNC Gene Name 1 | HGNC gene symbol associated with the first breakpoint of the control RNA fusion variant. |
| HGNC Gene Name 2 | HGNC gene symbol associated with the second breakpoint of the control RNA fusion variant. |

- **DNA NTC Library QC Metrics**—Contains information on the quality control metric that was evaluated for the DNA No-Template Control. The status of PASS indicates that the value for the metric is within the lower specification limit (LSL) and upper specification limit (USL) ranges. The status of FAIL indicates that value for the metric is outside of LSL or USL range. N/A values are listed if the DNA No-Template Control was not included in the sequencing run.

Table 20 DNA NTC Library QC Metrics

| Metric | Description | Units | Quality Threshold |
|----------------------|--|-------|-------------------|
| MEDIAN_EXON_COVERAGE | Median exon fragment coverage across all exon bases. | Count | ≤ 8 |

- **RNA NTC Library QC Metrics**—Contains information on the quality control metric that was evaluated for the RNA No-Template Control. The status of PASS indicates that the value for the metric is within the lower specification limit (LSL) and upper specification limit (USL) ranges. The status of FAIL indicates that value for the metric is outside of LSL or USL range. N/A values are listed if the RNA No-Template Control was not included in the sequencing run.

Table 21 RNA NTC Library QC Metrics

| Metric | Description | Units | Quality Threshold |
|--------------------------|--|-------|-------------------|
| GENE_ABOVE_MEDIAN_CUTOFF | The number of genes for which the median deduped read depth across all loci spanned for each gene is > 20. | Count | ≤ 1 |

Metrics Output

File name: `MetricsOutput.tsv`

The metrics output is a tab-delimited file that provides quality control information for patient samples that were included in the run.

The metrics output file contains the following sections and their associated fields:

- **Header**—Contains general information about the file and the run.

Table 22 Metrics Output File Header

| Field | Description |
|------------------|---|
| Output Date | Date this file was created. |
| Output Time | Time this file was created. |
| Workflow Version | The version of the analysis pipeline/workflow. |
| Module Version | The version of the TSO Comprehensive analysis module. |
| Run ID | The ID of the sequencing run. |
| Run Name | The name of the sequencing run. |

- **Run QC Metrics**—Contains quality control information for the sequencing run. This section corresponds to the Run QC status in the TSO Comprehensive (JP) report and contains one row per QC metric that contributes to Run QC status. All QC metrics in this section must pass for Run QC to pass. Refer to [Run Quality Control on page 8](#) for analysis details. Refer to [Quality Control Metrics on page 54](#) for metric descriptions and thresholds.

Table 23 Run QC Metrics

| Column | Description |
|--------------|---|
| Metric (UOM) | QC metric name and unit of measurement. |
| LSL | Lower specification limit (inclusive). |
| USL | Upper specification limit (inclusive). |
| Value | QC metric value. |
| PASS/FAIL | Indicates whether the sample passed or failed the quality control metric. Possible values include PASS, FAIL, or N/A. |

- **Analysis Status**—Contains information on whether analysis was completed for each patient sample, and whether any samples failed due to a software error. Each column in this section corresponds to a patient sample (Sample ID is used for the column name).

Table 24 Analysis Status

| Field | Description |
|---------------------|---|
| COMPLETED_ALL_STEPS | Indicates whether the sample completed all steps of the analysis. Possible values include TRUE and FALSE. If the value is FALSE, contact Illumina technical support for more information. |
| FAILED_STEPS | A list of any failed analysis steps due to a software error. Contact Illumina technical support for more information if any step is listed here. |
| STEPS_NOT_EXECUTED | A list of any analysis steps not executed due to a software error. Contact Illumina technical support for more information if any step is listed here. |

- **QC Metrics Sections for Patient Samples**—A section is included for each type of quality control used for patient samples. The following table notes where a quality control status in the TSO Comprehensive (JP) report corresponds to a section.

Table 25 QC Metrics Sections for Patient Samples

| Section | Description | Corresponding QC Category in TSO Comprehensive (JP) Report |
|--|---|--|
| DNA Library QC Metrics | QC metrics used as validity criteria for DNA sample libraries. Refer to Quality Control for DNA Sample Libraries on page 12 for analysis details. Refer to Quality Control Metrics on page 54 for metric descriptions and thresholds. | DNA Library QC |
| DNA Library QC Metrics for Small Variant Calling and TMB | QC metrics used as a validity criteria for small variants and TMB in a DNA sample library. Refer to Quality Control for DNA Sample Libraries on page 12 for analysis details. Refer to Quality Control Metrics on page 54 for metric descriptions and thresholds. | DNA Small Variant & TMB QC |

| Section | Description | Corresponding QC Category in TSO Comprehensive (JP) Report |
|--------------------------------|--|--|
| DNA Library QC Metrics for MSI | QC metrics used as validity criteria for MSI in a DNA Solid-FFPE sample library. Refer to Quality Control for DNA Sample Libraries on page 12 for analysis details. Refer to Quality Control Metrics on page 54 for metric descriptions and thresholds. | DNA MSI QC |
| DNA Library QC Metrics for CNV | QC metrics used as validity criteria for gene amplifications in a DNA Solid-FFPE sample library. Refer to Quality Control for DNA Sample Libraries on page 12 for analysis details. Refer to Quality Control Metrics on page 54 for metric descriptions and thresholds. | DNA Copy Number Variant QC |
| DNA Expanded Metrics | DNA Expanded Metrics are for information only and do not directly indicate the quality of DNA libraries. Refer to Quality Control for DNA Sample Libraries on page 12 for analysis details. Refer to DNA Expanded Metrics on page 58 for metric descriptions. | N/A |
| RNA Library QC Metrics | QC metrics used as validity criteria for RNA sample libraries. Refer to Quality Control for RNA Sample Libraries on page 15 for analysis details. Refer to Quality Control Metrics on page 54 for metric descriptions and thresholds. | RNA Library QC |
| RNA Expanded Metrics | RNA Expanded Metrics are for information only and do not directly indicate the quality of RNA libraries. Refer to Quality Control for RNA Sample Libraries on page 15 for analysis details. Refer to RNA Expanded Metrics on page 58 for metric descriptions and thresholds. | N/A |

Each section contains the following columns:

- Metric (UOM)—The QC metric name and unit of measurement.
- LSL—Lower specification limit (inclusive).
- USL—Upper specification limit (inclusive).
- One column per sample (named with Sample ID).

Each section contains the following rows:

- One row per QC metric.
- PASS/FAIL—Indicates whether the sample passed or failed for the type of quality control. A status of PASS indicates that the sample values for the metrics are within LSL and USL range. A status of FAIL indicates that sample values for one or more of the metrics are outside of the LSL or USL range. This row is not included for DNA Expanded Metrics or RNA Expanded Metrics.
- **Notes**—Contains a list of notes describing the content of the file.

Low Depth Report

File name: {SAMPLE_ID}_LowDepthReport.tsv

The low depth report is a tab-delimited file created for each patient sample. The file includes a listing of genomic position ranges with a total sequencing depth < 100 and for which a passing variant was not detected. These positions have insufficient sequencing depth to rule out the presence of a small variant. Positions on the block list are excluded from the report.

The low depth report is not regenerated during Report Regeneration.

The low depth report contains the following sections and their associated fields:

- **Header**—Contains general information about the file and the run.

Table 26 Header Information

| Field | Description |
|----------------------------------|---|
| Sample ID | Sample ID of the patient sample. |
| Report Date | The date the low depth report was generated. |
| Run ID | The ID of the sequencing run. |
| Run Date | The date of the sequencing run. |
| Knowledge base version | The version of the KB that was installed when the low depth report was generated. |
| Knowledge base published date | The date associated with KB that was installed when the low depth report was generated. |
| Local Run Manager Module version | The version of the TSO Comprehensive analysis module. |

- **Genomic Range List**—Contains a list of genomic position ranges with low depth. Contiguous genomic positions with low depth overlapping the same genes are combined into a single row.

Table 27 Genomic Range List

| Column | Description |
|--------|------------------------|
| Chrom | Chromosome. |
| Start | Start position (hg19). |

| Column | Description |
|--------|---|
| End | End position (hg19). |
| Gene | One or more gene symbols overlapping the genomic range based on the RefSeq database included in the KB. |

Output Folder Structure

This section describes the content of each output folder generated during analysis.

- IVD
 - IVD_Reports
 - {SampleID}_TSOCompJPModule_KB{version}_Report.pdf—TSO Comprehensive (JP) report (PDF format) per patient sample
 - {SampleID}_TSOCompJPModule_KB{version}_Report.json—TSO Comprehensive (JP) report (JSON format) per patient sample
 - {SampleID}_LowDepthReport.tsv—Low depth report per patient sample
 - MetricsOutput.tsv—Metrics output
 - ControlOutput.tsv—Control output report
- **Logs_Intermediates**—Logs and intermediate files generated during the analysis pipeline/workflow. Intermediate files are intended to help with troubleshooting only. The information contained in the intermediate files is not intended to be used for clinical reporting or patient management. Performance of any variants identified in these files, other than validated variants, has not been demonstrated. Validated variants are variants with demonstrated performance characteristics. Each folder represents one step of the analysis pipeline/workflow. The TSO Comprehensive analysis module appends RNA or DNA to the Sample ID folder names during processing.

View Analysis Results

1. From the Local Run Manager dashboard, select the run name.
2. From the Run Overview tab, review the sequencing run metrics.
3. To change the analysis data file location for future requeues of the selected run, select the **Edit** icon, and edit the output run folder file path.
The file path leading up to the output run folder is editable. The output run folder name cannot be changed.
4. **[Optional]** Select the **Copy to Clipboard** icon to copy the output run folder file path.
5. Select the Sequencing Information tab to review run parameters and consumables information.
6. Select the Samples & Results tab to view the analysis report.

- If analysis was requeued, select the appropriate analysis from the Select Analysis drop-down list.

7. [Optional] Select the **Copy to Clipboard** icon to copy the Analysis Folder file path.

Samples & Results

The Samples & Results screen displays the analysis results associated with the selected run and provides the option to reanalyze the run with different parameters. A table at the top of the screen provides the start date of the currently selected analysis run and the type of run (initial analysis, analysis requeue, or report regeneration).

Run Level Metrics

The Run Level Metrics section of the Samples & Results screen displays a run QC metric status of PASS or FAIL for each Run QC metric. Run QC metric statuses are sourced from the `MetricsOutput.tsv` file (refer to [Metrics Output on page 43](#)). Refer to [Quality Control Metrics on page 54](#) for metric descriptions and thresholds.

Controls

Controls are designated in the Run Setup screen of the TSO Comprehensive analysis module. Results for controls are displayed in the Controls section of the Samples & Results screen. The Controls section displays the following columns for each sample designated as a control:

- **Sample ID**
- **Type**—Control type. Possible values are DNA External Control, DNA No-Template Control, RNA External Control, and RNA No-Template Control. The installed KB does not affect the available control types.
- **Analysis Complete?**—Possible values are TRUE and FALSE. Controls marked as TRUE in the Analysis Complete? column have completed control analysis. If a control is marked FALSE, a software error has occurred. Contact Illumina technical support for more information.
- **Outcome**—Possible values are PASS and FAIL. DNA and RNA controls are evaluated independently. Refer to the following table for outcome value interpretation:

| Control type | Outcome | Interpretation |
|--------------|---------|--|
| DNA | PASS | Cross-contamination between libraries is not indicated. |
| No-Template | FAIL | Cross-contamination between libraries is indicated. DNA samples in the library preparation event and all associated sequencing runs are invalid. |

| Control type | Outcome | Interpretation |
|--------------------|---------|--|
| RNA No-Template | PASS | Cross-contamination between libraries is not indicated. |
| | FAIL | Cross-contamination between libraries is indicated. RNA samples in the library preparation event and all associated sequencing runs are invalid. |
| DNA External | PASS | Expected variants have been detected. |
| | FAIL | Variant calling specifications have not been met and DNA samples in the sequencing run are invalid. |
| RNA External | PASS | Expected variants have been detected. |
| | FAIL | Variant calling specifications have not been met and RNA samples in the sequencing run are invalid. |

Sample Level Metrics

The Sample Level Metrics section of the Samples & Results screen displays quality control information for patient samples that were included in the run. Patient sample quality control results are sourced from the `MetricsOutput.tsv` file (refer to [Metrics Output on page 43](#)). The Sample Level Metrics section displays the following columns for each patient sample:

- **Sample**—The sample ID.
- **Analysis Complete?**—Possible values are TRUE and FALSE. Samples marked as TRUE in the Analysis Complete? column have completed analysis successfully. If a sample is marked FALSE in this column, a software error has occurred. Contact Illumina technical support for more information.
- **DNA Library QC**—Possible values are PASS and FAIL. Indicates whether the sample passed or failed DNA library QC, which applies to the DNA library that was sequenced. Corresponds to DNA Library QC in the TSO Comprehensive (JP) report. A dash (–) is shown if a DNA library was not sequenced, or Run QC has a value of FAIL.
- **DNA Variants and Biomarkers**
 - **Small Variants and TMB**—Possible values are PASS and FAIL. Indicates whether the sample passed or failed QC for small variants and TMB in the DNA library. Corresponds to DNA Small Variant and TMB QC in the TSO Comprehensive (JP) report. A dash (–) is shown if a DNA library was not sequenced, Run QC has a value of FAIL, or DNA Library QC has a value of FAIL.
 - **MSI**—Possible value are PASS and FAIL. Indicates whether the sample passed or failed QC for MSI in the DNA library. Corresponds to DNA MSI QC in the TSO Comprehensive (JP) report. A dash (–) is shown if a DNA library was not sequenced, Run QC has a value of FAIL, or DNA Library QC has a value of FAIL.
 - **CNV**—Possible value are PASS and FAIL. Indicates whether the sample passed or failed QC for gene amplifications in the DNA Solid-FFPE library. Corresponds to DNA Copy Number Variant QC in the TSO Comprehensive (JP) report. A dash (–) is shown if a DNA Solid-FFPE library was not sequenced, Run QC has a value of FAIL, or DNA Library QC has a value of FAIL.

- **RNA Library QC**—Possible values are PASS and FAIL. Indicates whether the sample passed or failed RNA library QC, which applies to the RNA Solid-FFPE library that was sequenced. Corresponds to RNA Library QC in the TSO Comprehensive (JP) report. A dash (–) is shown if an RNA library was not sequenced, or Run QC has a value of FAIL.

Report Regeneration

Report regeneration allows one or more reports to be regenerated without repeating all secondary analysis steps.

Report regeneration is much faster than a full analysis requeue but has different features:

- **Scope**—Report regeneration rebuilds the TSO Comprehensive (JP) report but skips some analysis steps. You can change the sex or tumor type for one or more samples or install a new KB to produce a new report reflecting these changes. Each sample must be manually selected for report regeneration, while an analysis requeue automatically selects all samples by default. Individual samples can be removed for analysis requeue.
- **Analysis run status**—Report regeneration requires a successful analysis run as input, while analysis requeue can be used in scenarios where analysis has failed.
- **Editable fields**—Report regeneration allows changes to the Sex and Tumor Type fields, while analysis requeue allows any of the fields selected during run setup to be changed.
- **TSO Comprehensive analysis module version**—Report regeneration requires a successful analysis from a matching version of the TSO Comprehensive analysis module.
- **Workflow KB version**—Report regeneration requires a successful analysis using a KB with a matching version of the RefSeq database.
- **Run Input Settings**—Report regeneration run inputs are automatically set to the values from the most recent successful secondary analysis run. The run inputs for an analysis requeue are automatically set to the values from the most recent analysis attempt (including failed analysis runs).

This feature is only accessible to TSO Comprehensive analysis module admin users or a non-admin user with requeue analysis permissions assigned. For more information on TSO Comprehensive analysis module user management, refer to *NextSeq 550Dx Instrument Reference Guide for Japan (document # 1000000009513)*.

Regenerate a Report or Requeue Analysis

1. From the run dashboard, locate a run with a status of Analysis Completed. Select the vertical ellipses icon and select **Requeue**.

Relinking runs that have been deleted from the local temp folder is required to requeue analysis. For more information on TSO Comprehensive analysis module user management, refer to *NextSeq 550Dx Instrument Reference Guide for Japan (document # 1000000009513)*.

2. Select **Edit Setup** in the Requeue Analysis pop-up.
3. Use the dropdown at the top of the Requeue Analysis screen to select report regeneration or full analysis requeue.

NOTE Always review run inputs for each sample before saving a run. Report regeneration run inputs are automatically set to the values from the most recent successful secondary analysis run.

4. Samples from the previously completed run are displayed in a table. Use the + buttons on the right of the table to mark desired samples for report regeneration. All samples in a run are excluded from report regeneration by default and must be added individually. Report regeneration is not available for samples originally analyzed as controls, which require full analysis requeue.
5. When all desired samples have been marked for report regeneration, select **Requeue Analysis**.

Viewing Report Regeneration Results

Regenerated reports for samples marked for report regeneration can be viewed along with other completed analyses in the Samples and Runs screen in TSO Comprehensive analysis module. Reports produced using report regeneration are marked as Report Regeneration in the Analysis Type field at the top of the Samples and Runs screen.

Appendix A QC Metrics Flowchart


The following flowchart describes the QC metrics that are listed on the TSO Comprehensive (JP) report. If Run QC fails, then no other QC steps are assessed, and all are marked as N/A. If DNA or RNA is not sequenced or fail Library QC, then any corresponding variant types are not included in Tumor Profiling results. DNA Library QC is a measure of contamination. If it does not pass, then the downstream DNA QC Metrics (DNA MSI QC, DNA small variants & TMB QC, and DNA CNV QC) are marked as N/A. For more information, refer to the following sections and tables:

- [Analysis Methods on page 7](#)
- [Quality Control on page 20](#)
- [Run QC Metrics on page 43](#)
- [Quality Control for DNA Sample Libraries on page 12](#)
- [Sample Level Metrics on page 49](#)
- [Appendix B QC Metrics on page 54](#)

The flowchart does not map the controls. The results from the controls do not impact the QC metrics on the TSO Comprehensive (JP) PDF or JSON report. Failure of controls invalidates sample results separate of QC results as described in [TruSight Oncology Comprehensive \(JP\) Report on page 18](#). The use of controls is described in [Controls on page 5](#). For additional controls information, refer to the *TruSight Oncology Comprehensive (JP) Assay Workflow Guide (document # 200041566)*.

The flowchart does not map the position-level QC results. Position-level QC results for the Tumor Profiling section are provided in the Low Depth Report (refer to [Low Depth Reporting for DNA Sample Libraries on page 12](#)).

Figure 2 QC Metrics Flowchart



Appendix B QC Metrics

Quality Control Metrics

Table 28 TSO Comprehensive Report Result QC Metrics

| Output Type | Metric | Specification | Description | Impact of Specification Failure* |
|----------------|------------------|---------------|---|--|
| Sequencing Run | PCT_PF_READS (%) | ≥ 80.0 | Percentage of reads passing filter (PF). | Sequencing run invalidated, no results reported for any sample in the run. |
| | PCT_Q30_R1 (%) | ≥ 80.0 | Average percent of base calls with quality score of Q30 or higher for Read 1. | |
| | PCT_Q30_R2 (%) | ≥ 80.0 | Average percent of base calls with quality score of Q30 or higher for Read 2. | |

| Output Type | Metric | Specification | Description | Impact of Specification Failure* |
|---------------|---------------------|--|--|----------------------------------|
| DNA Libraries | CONTAMINATION_SCORE | ≤ 3106 OR > 3106 and P_VALUE ≤ 0.049 | A metric assessing the likelihood of contamination using the VAF of common variants. The contamination score is based on VAF distribution of SNPs. The contamination P value used to assess highly rearranged genomes, only applicable when contamination score is above Upper Spec Limit. | No DNA results reported. |

| Output Type | Metric | Specification | Description | Impact of Specification Failure* |
|-------------|-------------------------------------|---------------|--|---|
| | MEDIAN_INSERT_SIZE (bp) | ≥ 70 | The median fragment length in the sample. | No TMB or small DNA variant results reported. |
| | MEDIAN_EXON_COVERAGE (count) | ≥ 150 | Median exon fragment coverage across all exon bases. | |
| | PCT_EXON_50X (%) | ≥ 90.0 | Percent exon bases with 50X fragment coverage. | |
| | USABLE_MSI_SITES (count) | ≥ 40 | The number of MSI sites usable for MSI calling (Number of microsatellite sites with sufficient spanning reads to identify microsatellite instability). | No MSI results reported. |
| | COVERAGE_MAD (count) | ≤ 0.210 | The median of absolute deviations from the median of the normalized count of each CNV target region. | No gene amplification results reported. |
| | MEDIAN_BIN_COUNT_CNV_TARGET (count) | ≥ 1.0 | The median raw bin count per CNV target. | |

| Output Type | Metric | Specification | Description | Impact of Specification Failure* |
|---------------|-----------------------------------|------------------|--|---|
| RNA Libraries | MEDIAN_INSERT_SIZE (bp) | ≥ 80 | The median fragment length in the sample. | No fusion or splice variant results reported. |
| | MEDIAN_CV_GENE_500X (coefficient) | ≤ 0.93 | MEDIAN_CV_GENE_500X is a measure of coverage uniformity. For each gene with at least 500x coverage, the coefficient of variation in coverage across the gene body is computed. This metric is the median of these values. A high value indicates a high level of variation and indicates a problem in library preparation such as low sample input and/or probe pulldown issues. This metric is computed using all reads (including reads marked as duplicates). | |
| | TOTAL_ON_TARGET_READS (count) | $\geq 9,000,000$ | The total number of reads that map to the target regions. This metric is computed using all reads (including reads marked as duplicates). | |

*Successful results show PASS.

DNA Expanded Metrics

DNA expanded metrics are provided for information only. They can be informative for troubleshooting but are provided without explicit specification limits and are not directly used for sample quality control. For additional guidance, contact Illumina Technical Support.


| Metric | Description | Units |
|------------------------|--|-------|
| TOTAL_PF_READS | Total reads passing filter. | Count |
| MEAN_FAMILY_SIZE | The sum of the reads in each family divided by the number of families after correction, collapsing, and filtering on supporting reads. | Count |
| MEDIAN_TARGET_COVERAGE | The median coverage of bases. | Count |
| PCT_CHIMERIC_READS | Percent of chimeric reads. | % |
| PCT_EXON_100X | Percent of exon bases with greater than 100X coverage. | % |
| PCT_READ_ENRICHMENT | Percentage of reads that cross any part of the target region vs total reads. | % |
| PCT_USABLE_UMI_READS | The percentage of reads with usable UMIs. | % |
| MEAN_TARGET_COVERAGE | The mean coverage of bases. | Count |
| PCT_ALIGNED_READS | Percent of reads that aligned to the reference genome. | % |
| PCT_CONTAMINATION_EST | Percent of contamination of the sample. | % |
| PCT_PF_UQ_READS | Percent unique reads passing filter. | % |
| PCT_TARGET_0.4X_MEAN | Percent target bases with target coverage greater than 0.4 times the mean. | % |
| PCT_TARGET_100X | Percent target bases with greater than 100X coverage. | % |
| PCT_TARGET_250X | Percent target bases with greater than 250X coverage. | % |

RNA Expanded Metrics

RNA expanded metrics are provided for information only. They can be informative for troubleshooting but are provided without explicit specification limits and are not directly used for sample quality control. For additional guidance, contact Illumina Technical Support.

| Metric | Description | Units |
|-------------------------------------|---|--------------|
| PCT_ CHIMERIC_ READS | Percentage of reads that are aligned as two segments that map to non-consecutive regions in the genome. | % |
| PCT_ON_ TARGET_ READS | Percentage of reads that cross any part of the target region vs total reads. A read that partially maps to a target region is counted as on target. | % |
| SCALED_ MEDIAN_GENE_ COVERAGE | Median of median base coverage of genes scaled by length. An indication of median coverage depth of genes in the panel. | Count |
| TOTAL_PF_ READS | Total number of reads passing filter. | Count |

Appendix C TSO Comprehensive (JP) Report Reference



- A. Refer to [Appendix A QC Metrics Flowchart on page 52](#) for details.
- B. The Other Alterations and Biomarkers Identified section contains tumor profiling information. Associations can be due to therapeutic, diagnostic, or prognostic evidence.
- C. According to the KB, there is evidence of clinical significance for this biomarker in this cancer type based on information from approved therapies, clinical guidelines, or both. For more information, refer to [Genomic Findings with Evidence of Clinical Significance on page 16](#) and the Level 2 table [Genomic Findings with Evidence of Clinical Significance on page 25](#).
- D. According to the KB, there is limited or no clinical evidence for a genomic finding within the cancer type. There might be preclinical data or data in other cancer types where the biomarker is predictive of response to an approved or investigational therapy. For more information, refer to [Genomic Findings with Potential Clinical Significance on page 16](#) and the Level 3 table [Genomic Findings with Potential Clinical Significance on page 27](#).
- E. TMB and MSI are listed in Level 3 Genomic Findings with Potential Clinical Significance. Refer to [Tumor Mutational Burden on page 11](#) and [Microsatellite Instability Status on page 12](#).
- F. If there are two variants listed in a single row (not pictured), there is clinical meaning for these variants when they are detected together. Resistance mutations or other sources can be the cause. Refer to examples in [Tumor Profiling of Variants on page 15](#).

Appendix D MNVs, Indels, and Deletions in EGFR and RET Detectable by Phased Variant Caller

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|-----------------------|------------------|------|---|
| chr7 | 55242462 | CAAGGAATTAAGAGAA | C | EGFR | NP_005219.2:p.(Lys745_Glu749del) |
| chr7 | 55242463 | AAGGAATTAAGAGAAG | A | EGFR | NP_005219.2:p.(Lys745_Ala750delinsThr) |
| chr7 | 55242464 | AGGAATTAAGAGA | A | EGFR | NP_005219.2:p.(Glu746_Glu749del) |
| chr7 | 55242464 | AGGAATTAAGAGAAGC | A | EGFR | NP_005219.2:p.(Glu746_Ala750del) |
| chr7 | 55242465 | GGAATTAAGA | G | EGFR | NP_005219.2:p.(Leu747_Glu749del) |
| chr7 | 55242465 | GGAATTAAGAGAAG | AATTC | EGFR | NP_005219.2:p.(Glu746_Ala750delinsIlePro) |
| chr7 | 55242465 | GGAATTAAGAGAAGCAA | AATTC | EGFR | NP_005219.2:p.(Glu746_Thr751delinsIlePro) |
| chr7 | 55242465 | GGAATTAAGAGAAGCAAC | AAT | EGFR | NP_005219.2:p.(Glu746_Thr751delinsIle) |
| chr7 | 55242465 | GGAATTAAGAGAAGCAACA | G | EGFR | NP_005219.2:p.(Glu746_Thr751del) |
| chr7 | 55242465 | GGAATTAAGAGAAGCAACATC | AAT | EGFR | NP_005219.2:p.(Glu746_Ser752delinsIle) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|-----------------------|------------------|------|---|
| chr7 | 55242465 | GGAATTAAGAGAAGCA | G | EGFR | NP_005219.2:p.(Glu746_Ala750del) |
| chr7 | 55242466 | GAATTAAGAGAAGCAACAT | G | EGFR | NP_005219.2:p.(Glu746_Ser752delinsAla) |
| chr7 | 55242466 | GAATTAAGAGAAGCAA | G | EGFR | NP_005219.2:p.(Glu746_Thr751delinsAla) |
| chr7 | 55242467 | AATTAAGAGAAGCAAC | A | EGFR | NP_005219.2:p.(Leu747_Thr751del) |
| chr7 | 55242467 | AATTAAGAGAAGCAACATC | A | EGFR | NP_005219.2:p.(Glu746_Ser752delinsAsp) |
| chr7 | 55242467 | AATTAAGAGAAGCAACATC | T | EGFR | NP_005219.2:p.(Glu746_Ser752delinsVal) |
| chr7 | 55242467 | AATTAAGAGAAGCAACATCTC | TCT | EGFR | NP_005219.2:p.(Glu746_Pro753delinsValSer) |
| chr7 | 55242467 | AATTAAGAGAAGCAACA | TTGCT | EGFR | NP_005219.2:p.(Glu746_Thr751delinsValAla) |
| chr7 | 55242467 | AATTAAGAGAAGCAAC | T | EGFR | NP_005219.2:p.(Glu746_Thr751delinsVal) |
| chr7 | 55242468 | ATTAAGAGAAGCAACATCT | A | EGFR | NP_005219.2:p.(Leu747_Ser752del) |
| chr7 | 55242468 | ATTAAGAGAAGCAAC | GCA | EGFR | NP_005219.2:p.(Leu747_Thr751delinsGln) |
| chr7 | 55242468 | ATTAAGAGAAG | GC | EGFR | NP_005219.2:p.(Leu747_Ala750delinsPro) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|---------------------------|------------------|------|--|
| chr7 | 55242469 | TTAAGAGAAG | C | EGFR | NP_005219.2:p.(Leu747_ Ala750delinsPro) |
| chr7 | 55242469 | TTAAGAGAAGCAA | C | EGFR | NP_005219.2:p.(Leu747_ Thr751delinsPro) |
| chr7 | 55242469 | TTAAGAGAAGCAACATCT | CAA | EGFR | NP_005219.2:p.(Leu747_ Ser752delinsGln) |
| chr7 | 55242469 | TTAAGAGAAGCAACATCTCC | CA | EGFR | NP_005219.2:p.(Leu747_ Pro753delinsGln) |
| chr7 | 55242469 | TTAAGAGAAGCAACATCTC | T | EGFR | NP_005219.2:p.(Leu747_ Pro753delinsSer) |
| chr7 | 55242469 | TTAAGAGAAGCAA | T | EGFR | NP_005219.2:p.(Leu747_ Thr751delinsSer) |
| chr7 | 55242482 | CATCTCCGAAAGCCAACAAGGAAAT | C | EGFR | NP_005219.2:p.(Ser752_ Ile759del) |
| chr7 | 55249011 | AC | CCAGCGTGGAT | EGFR | NP_005219.2:p.(Ala767_ Val769dup) |
| chr10 | 43604549 | CTCAGACTTCCAGGGCCCAGGA | G | RET | NP_066124.1:p.(Asp378_ Gly385delinsGlu) |
| chr10 | 43609928 | ATCCTACTGTGCGACGAGCTG | CACAC | RET | NP_066124.1:p.(Asp627_ Leu633delinsAlaHis) |
| chr10 | 43609928 | ATCCTACTGTGCGACGAGCTG | CACAT | RET | NP_066124.1:p.(Asp627_ Leu633delinsAlaHis) |
| chr10 | 43609928 | ATCCTACTGTGCGACGAGCTG | CCCAC | RET | NP_066124.1:p.(Asp627_ Leu633delinsAlaHis) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|--------------------------------|------------------|------|---|
| chr10 | 43609928 | ATCCACTGTGCGACGAGCTG | CCCAT | RET | NP_066124.1:p.(Asp627_Leu633delinsAlaHis) |
| chr10 | 43609928 | ATCCACTGTGCGACGAGCTG | CGCAC | RET | NP_066124.1:p.(Asp627_Leu633delinsAlaHis) |
| chr10 | 43609928 | ATCCACTGTGCGACGAGCTG | CGCAT | RET | NP_066124.1:p.(Asp627_Leu633delinsAlaHis) |
| chr10 | 43609928 | ATCCACTGTGCGACGAGCTG | CTCAC | RET | NP_066124.1:p.(Asp627_Leu633delinsAlaHis) |
| chr10 | 43609928 | ATCCACTGTGCGACGAGCTG | CTCAT | RET | NP_066124.1:p.(Asp627_Leu633delinsAlaHis) |
| chr10 | 43609933 | CTGTGCGACGAGCTGTGCCGCACGGTGATC | TGCGAT | RET | NP_066124.1:p.(Leu629_Ile638delinsCysAsp) |
| chr10 | 43609933 | CTGTGCGACGAGCTGTGCCGCACGGTGATC | TGTGAT | RET | NP_066124.1:p.(Leu629_Ile638delinsCysAsp) |
| chr10 | 43609933 | CTGTGCGACGAGCTGTGCCGCACGGTGAT | TGCGA | RET | NP_066124.1:p.(Leu629_Ile638delinsCysAsp) |
| chr10 | 43609933 | CTGTGCGACGAGCTGTGCCGCACGGTGAT | TGTGA | RET | NP_066124.1:p.(Leu629_Ile638delinsCysAsp) |
| chr10 | 43609936 | TGC | GCT | RET | NP_066124.1:p.(Cys630Ala) |
| chr10 | 43609940 | ACGAGCTG | TA | RET | NP_066124.1:p.(Asp631_Leu633delinsVal) |
| chr10 | 43609940 | ACGAGCTG | TC | RET | NP_066124.1:p.(Asp631_Leu633delinsVal) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|-------------------------|------------------|------|---|
| chr10 | 43609940 | ACGAGCTGTGCCGCACGGTGAT | C | RET | NP_066124.1:p.(Asp631_Ile638delinsAla) |
| chr10 | 43609940 | ACGAGCTGTGCCGCACGGTGATC | CA | RET | NP_066124.1:p.(Asp631_Ile638delinsAla) |
| chr10 | 43609940 | ACGAGCTGTGCCGCACGGTGATC | CG | RET | NP_066124.1:p.(Asp631_Ile638delinsAla) |
| chr10 | 43609940 | ACGAGCTGTGCCGCACGGTGATC | CT | RET | NP_066124.1:p.(Asp631_Ile638delinsAla) |
| chr10 | 43609940 | ACGAGCTG | TT | RET | NP_066124.1:p.(Asp631_Leu633delinsVal) |
| chr10 | 43609941 | CGAGCTG | A | RET | NP_066124.1:p.(Asp631_Leu633delinsGlu) |
| chr10 | 43609942 | GAGCTGTGCCGCA | AGCT | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCA | AGTT | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | AGCAGC | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | AGCAGT | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | AGCTCA | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | AGCTCC | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|--------------------------|------------------|------|---|
| chr10 | 43609942 | GAGCTGTGCCGCACG | AGCTCT | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | AGTAGC | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | AGTAGT | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | AGTTCA | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | AGTTCC | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | AGTTCT | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACGGTGATCGCA | CACAGG | RET | NP_066124.1:p.(Glu632_Ala639delinsHisArg) |
| chr10 | 43609942 | GAGCTGTGCCGCACGGTGATCGCA | CACCGC | RET | NP_066124.1:p.(Glu632_Ala639delinsHisArg) |
| chr10 | 43609942 | GAGCTGTGCCGCACGGTGATCGCA | CACCGG | RET | NP_066124.1:p.(Glu632_Ala639delinsHisArg) |
| chr10 | 43609942 | GAGCTGTGCCGCACGGTGATCGCA | CACCGT | RET | NP_066124.1:p.(Glu632_Ala639delinsHisArg) |
| chr10 | 43609942 | GAGCTGTGCCGCACGGTGATCGCA | CATAGG | RET | NP_066124.1:p.(Glu632_Ala639delinsHisArg) |
| chr10 | 43609942 | GAGCTGTGCCGCACGGTGATCGCA | CATCGC | RET | NP_066124.1:p.(Glu632_Ala639delinsHisArg) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|--------------------------|------------------|------|---|
| chr10 | 43609942 | GAGCTGTGCCGCACGGTGATCGCA | CATCGG | RET | NP_066124.1:p.(Glu632_Ala639delinsHisArg) |
| chr10 | 43609942 | GAGCTGTGCCGCACGGTGATCGCA | CATCGT | RET | NP_066124.1:p.(Glu632_Ala639delinsHisArg) |
| chr10 | 43609942 | GAGCTGTGCCGCACGGTGATCGC | CACAG | RET | NP_066124.1:p.(Glu632_Ala639delinsHisArg) |
| chr10 | 43609942 | GAGCTGTGCCGCACGGTGATCGC | CACCG | RET | NP_066124.1:p.(Glu632_Ala639delinsHisArg) |
| chr10 | 43609942 | GAGCTGTGCCGCACGGTGATCGC | CATAG | RET | NP_066124.1:p.(Glu632_Ala639delinsHisArg) |
| chr10 | 43609942 | GAGCTGTGCCGCACGGTGATCGC | CATCG | RET | NP_066124.1:p.(Glu632_Ala639delinsHisArg) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | TCAAGC | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | TCAAGT | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | TCATCA | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | TCATCC | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | TCATCT | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | TCCAGC | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|------------------|------------------|------|---|
| chr10 | 43609942 | GAGCTGTGCCGCACG | TCCAGT | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | TCCTCA | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | TCCTCC | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | TCCTCT | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | TCGAGC | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | TCGAGT | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | TCGTCA | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | TCGTCC | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | TCGTCT | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | TCTAGC | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | TCTAGT | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | TCTTCA | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|------------------|------------------|------|---|
| chr10 | 43609942 | GAGCTGTGCCGCACG | TCTTCC | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCACG | TCTTCT | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCA | TCAT | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCA | TCCT | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCA | TCGT | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609942 | GAGCTGTGCCGCA | TCTT | RET | NP_066124.1:p.(Glu632_Thr636delinsSerSer) |
| chr10 | 43609943 | AGCTG | TA | RET | NP_066124.1:p.(Glu632_Leu633delinsVal) |
| chr10 | 43609943 | AGCTG | TC | RET | NP_066124.1:p.(Glu632_Leu633delinsVal) |
| chr10 | 43609943 | AGCTGTGCCGCACGGT | CAGC | RET | NP_066124.1:p.(Glu632_Val637delinsAlaAla) |
| chr10 | 43609943 | AGCTGTGCCGCACGGT | CCGC | RET | NP_066124.1:p.(Glu632_Val637delinsAlaAla) |
| chr10 | 43609943 | AGCTGTGCCGCACGGT | CGGC | RET | NP_066124.1:p.(Glu632_Val637delinsAlaAla) |
| chr10 | 43609943 | AGCTGTGCCGCACGGT | CTGC | RET | NP_066124.1:p.(Glu632_Val637delinsAlaAla) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|----------------------------|------------------|------|--|
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TAAGACCA | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TAAGACCG | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TAAGACCT | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TAAGGCCA | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TAAGGCCT | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TACGACCA | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TACGACCG | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TACGACCT | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TACGCCCA | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TACGCCCG | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TACGCCCT | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TACGGCCA | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|----------------------------|------------------|------|--|
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TACGGCCG | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TACGGCCT | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TACGTCCA | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TACGTCCG | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TACGTCCT | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TCAGACCA | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TCAGACCG | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TCAGACCT | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TCAGGCCA | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TCAGGCCG | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TCAGGCCT | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TCCGACCA | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|----------------------------|------------------|------|--|
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TCCGACCG | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TCCGACCT | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TCCGCCCA | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TCCGCCCG | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TCCGCCCT | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TCCGGCCA | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TCCGGCCG | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TCCGGCCT | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TCCGTCCA | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TCCGTCCG | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TCCGTCTT | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TGAGACCA | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|----------------------------|------------------|------|--|
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TGAGACCG | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TGAGACCT | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TGAGGCCA | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TGAGGCCG | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TGAGGCCT | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TGCGACCA | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TGCGACCG | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TGCGACCT | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TGCGCCCA | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TGCGCCCG | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TGCGCCCT | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TGCGGCCA | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|----------------------------|------------------|------|--|
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TGCGGCCT | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TGCGTCCA | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TGCGTCCG | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TGCGTCCT | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TTAGACCA | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TTAGACCG | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TTAGACCT | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TTAGGCCA | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TTAGGCCG | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TTAGGCCT | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TTCGACCA | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TTCGACCG | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|----------------------------|------------------|------|--|
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TTCGACCT | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TTCGCCCA | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TTCGCCCG | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TTCGCCCT | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TTCGGCCA | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TTCGGCCG | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TTCGGCCT | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TTCGTCCA | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TTCGTCCG | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAGCC | TTCGTCTT | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAG | TAAGAC | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAG | TAAGGC | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|--------------------------|------------------|------|--|
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAG | TACGAC | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAG | TACGCC | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAG | TACGGC | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAG | TACGTC | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAG | TCAGAC | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAG | TCAGGC | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAG | TCCGAC | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAG | TCCGCC | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAG | TCCGGC | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAG | TCCGTC | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAG | TGAGAC | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAG | TGAGGC | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|--------------------------|------------------|------|--|
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAG | TGCGAC | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAG | TGCGCC | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAG | TGCGGC | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAG | TGCGTC | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAG | TTAGAC | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAG | TTAGGC | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAG | TTCGAC | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAG | TTCGCC | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAG | TTCGGC | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTGATCGCAG | TTCGTC | RET | NP_066124.1:p.(Glu632_Ala640delinsValArgPro) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTG | CAGCA | RET | NP_066124.1:p.(Glu632_Val637delinsAlaAla) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTG | CAGCC | RET | NP_066124.1:p.(Glu632_Val637delinsAlaAla) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|-------------------|------------------|------|--|
| chr10 | 43609943 | AGCTGTGCCGCACGGTG | CAGCT | RET | NP_066124.1:p.(Glu632_Val637delinsAlaAla) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTG | CCGCA | RET | NP_066124.1:p.(Glu632_Val637delinsAlaAla) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTG | CCGCC | RET | NP_066124.1:p.(Glu632_Val637delinsAlaAla) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTG | CCGCT | RET | NP_066124.1:p.(Glu632_Val637delinsAlaAla) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTG | CGGCA | RET | NP_066124.1:p.(Glu632_Val637delinsAlaAla) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTG | CGGCC | RET | NP_066124.1:p.(Glu632_Val637delinsAlaAla) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTG | CGGCT | RET | NP_066124.1:p.(Glu632_Val637delinsAlaAla) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTG | CTGCA | RET | NP_066124.1:p.(Glu632_Val637delinsAlaAla) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTG | CTGCC | RET | NP_066124.1:p.(Glu632_Val637delinsAlaAla) |
| chr10 | 43609943 | AGCTGTGCCGCACGGTG | CTGCT | RET | NP_066124.1:p.(Glu632_Val637delinsAlaAla) |
| chr10 | 43609943 | AGCTG | TT | RET | NP_066124.1:p.(Glu632_Leu633delinsVal) |
| chr10 | 43609944 | GCTGT | CGTAC | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|------------------|------------------|------|--|
| chr10 | 43609944 | GCTGT | CGTCC | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGT | CGTGC | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGT | CGTTC | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | CGTAAGA | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | CGTAAGG | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | CGTACGA | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | CGTACGG | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | CGTACGT | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | CGTCAGA | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | CGTCAGG | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | CGTCCGA | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | CGTCCGG | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|------------------|------------------|------|--|
| chr10 | 43609944 | GCTGTGC | CGTCCGT | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | CGTGAGA | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | CGTGAGG | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | CGTGCGA | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | CGTGCGG | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | CGTGCGT | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | CGTTAGA | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | CGTTAGG | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | CGTTCGA | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | CGTTCGG | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | CGTTCGT | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | TGTAAGA | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|------------------|------------------|------|--|
| chr10 | 43609944 | GCTGTGC | TGTAAGG | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | TGTACGA | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | TGTACGG | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | TGTACGT | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | TGTCAGA | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | TGTCCGA | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | TGTCCGG | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | TGTCCGT | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | TGTGAGA | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | TGTGAGG | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | TGTGCGA | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | TGTGCGG | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|------------------|------------------|------|--|
| chr10 | 43609944 | GCTGTGC | TGTGCGT | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | TGTTAGA | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | TGTTAGG | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | TGTTCGA | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | TGTTCGG | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGTGC | TGTTCGT | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGT | TGTAC | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGT | TGTCC | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGT | TGTGC | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609944 | GCTGT | TGTTC | RET | NP_066124.1:p.(Glu632_Cys634delinsAspValArg) |
| chr10 | 43609945 | CTGTGC | GTATGG | RET | NP_066124.1:p.(Leu633_Cys634delinsValTrp) |
| chr10 | 43609945 | CTGTGC | GTCTGG | RET | NP_066124.1:p.(Leu633_Cys634delinsValTrp) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|------------------|------------------|------|---|
| chr10 | 43609945 | CTGTGC | GTGTGG | RET | NP_066124.1:p.(Leu633_Cys634delinsValTrp) |
| chr10 | 43609945 | CTGTGC | GTTTGG | RET | NP_066124.1:p.(Leu633_Cys634delinsValTrp) |
| chr10 | 43609948 | TGC | CCA | RET | NP_066124.1:p. (Cys634Pro) |
| chr10 | 43609948 | TGC | CCG | RET | NP_066124.1:p. (Cys634Pro) |
| chr10 | 43609950 | CCGC | GGGA | RET | NP_066124.1:p.(Cys634_Arg635delinsTrpGly) |
| chr10 | 43609950 | CCGC | GGGG | RET | NP_066124.1:p.(Cys634_Arg635delinsTrpGly) |
| chr10 | 43609950 | CCGC | GGGT | RET | NP_066124.1:p.(Cys634_Arg635delinsTrpGly) |
| chr10 | 43609950 | CCGC | TCCAAAAGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCAAAAGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCAAAACGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCAAAACGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCAAAACGT | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|------------------|------------------|------|--|
| chr10 | 43609950 | CCGC | TCCAAAGAGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCAAAGAGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCAAAGCGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCAAAGCGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCAAAGCGT | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCCAAAGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCCAAAGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCCAAACGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCCAAACGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCCAAACGT | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCAAGAGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCAAGAGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|------------------|------------------|------|--|
| chr10 | 43609950 | CCGC | TCCAAGCGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCAAGCGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCAAGCGT | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCGAAAAGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCGAAAAGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCGAAACGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCGAAACGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCGAAACGT | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCGAAGAGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCGAAGAGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCGAAGCGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCGAAGCGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|------------------|------------------|------|--|
| chr10 | 43609950 | CCGC | TCCGAAGCGT | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCTAAAAGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCTAAAAGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCTAAACGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCTAAACGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCTAAACGT | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCTAAGAGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCTAAGAGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCTAAGCGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCTAAGCGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | CCGC | TCCTAAGCGT | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | C | TCCAAAA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|------------------|------------------|------|--|
| chr10 | 43609950 | C | TCCAAAG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | C | TCCCAA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | C | TCCCAAG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | C | TCCGAAA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | C | TCCGAAG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | C | TCCTAAA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609950 | C | TCCTAAG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CAAAAAGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CAAAACGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CAAAACGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CAAAACGT | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CAAAGAGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|------------------|------------------|------|--|
| chr10 | 43609952 | GC | CAAAGCGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CAAAGCGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CAAAGCGT | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CCAAAAGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CCAAACGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CCAAACGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CCAAACGT | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CCAAGAGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CCAAGAGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CCAAGCGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CCAAGCGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CCAAGCGT | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|------------------|------------------|------|--|
| chr10 | 43609952 | GC | CGAAAAGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CGAAAAGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CGAAACGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CGAAACGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CGAAACGT | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CGAAGAGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CGAAGAGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CGAAGCGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CGAAGCGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CGAAGCGT | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CTAAAAGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CTAAACGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|------------------|------------------|------|---|
| chr10 | 43609952 | GC | CTAAACGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CTAAACGT | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CTAAGAGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CTAAGCGA | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CTAAGCGG | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43609952 | GC | CTAAGCGT | RET | NP_066124.1:p.(Cys634_Arg635insProLys) |
| chr10 | 43613904 | TTG | ACT | RET | NP_066124.1:p. (Leu790Thr) |
| chr10 | 43615630 | TTCC | ACCA | RET | NP_066124.1:p.(Asp903_Ser904delinsGluPro) |
| chr10 | 43615630 | TTCC | ACCG | RET | NP_066124.1:p.(Asp903_Ser904delinsGluPro) |
| chr10 | 43615630 | TTCC | ACCT | RET | NP_066124.1:p.(Asp903_Ser904delinsGluPro) |
| chr10 | 43615630 | TTCC | GCCA | RET | NP_066124.1:p.(Asp903_Ser904delinsGluPro) |
| chr10 | 43615630 | TTCC | GCCG | RET | NP_066124.1:p.(Asp903_Ser904delinsGluPro) |

| Chromosome | Position (hg19) | Reference Allele | Alternate Allele | Gene | Amino Acid Change |
|------------|-----------------|------------------|------------------|------|---|
| chr10 | 43615630 | TTCC | GCCT | RET | NP_066124.1:p.(Asp903_Ser904delinsGluPro) |

Appendix E Install a Knowledge Base

The TSO Comprehensive analysis module requires an installed Knowledge Base (KB) to perform analysis. Illumina periodically releases new KBs. A Field Applications Scientist provides the KB and related files to the customer using a data sharing folder. When updating a KB, the previously installed KB is removed during the installation process. Do not install a KB while a sequencing run, analysis, or other installation process is in progress.



CAUTION

To avoid data loss, make sure that no other processes are in progress before following the installation instructions.



CAUTION

Navigating away from the Modules & Manifests page or closing the browser during the KB installation cancels the installation process.

1. Download the desired KB (*.zip) to a local directory on your instrument or a networked computer. Drive D is the preferred location.
2. Perform KB checksum verification as follows.
 - a. Perform a Windows search for PowerShell. Right-click on the program and select **Run as Administrator**.
 - b. Enter `Get-FileHash <KB file path>\<kbfilename.zip> -Algorithm MD5` in a PowerShell window to generate the MD5 checksum for the KB (for example, `Get-FileHash C:\Users\jdoe\Downloads\KB_EU_8.25.0.0325.zip -Algorithm MD5`).
 - c. Compare the output MD5 checksum against the KB checksum from the data sharing folder (for example, MD5: c2948cffffa0e1a891630f1ad5504046). If the checksums do not match, delete this KB file and download it again from the data sharing folder.
3. Open TSO Comprehensive analysis module on your instrument or the networked computer (local area network).
4. Sign in as an admin or a non-admin user with permission to edit module settings. For more information on TSO Comprehensive analysis module user management, refer to *NextSeq 550Dx Instrument Reference Guide for Japan (document # 1000000009513)*.
5. Use the Tools menu to navigate to the Modules & Manifests screen.
6. Select **TSO Comp (JP)**.
7. Select **Install New** under the Knowledge Base Version section. An installation wizard prompts you to browse to the location of the KB ZIP file.
8. Select the KB that was downloaded in step 1. The wizard also displays information about the KB including the name, version, RefSeq database version, and published date.

9. Select **Continue** in the installation wizard.

The installer verifies that the KB is compatible with the TSO Comprehensive analysis module and that the KB is not corrupt. It is not possible to launch a new TSO Comprehensive (JP) analysis while installing the KB. After installation is complete, the new KB is listed on the Modules & Manifests screen. The KB name and version are also displayed on the Create Run, Requeue Analysis, and Edit Run screens.

Appendix F Cybersecurity

Antivirus or Antimalware Software

The following antivirus (AV) or antimalware (AM) software products are compatible with NOS and TSO Comprehensive analysis module when configured following the instructions provided in the *NextSeq 550Dx Instrument Site Prep Guide for Japan (document # 1000000058597)*:

- Windows Defender/Windows Security
- BitDefender
- CrowdStrike

For additional details regarding network, firewall, and storage configurations, contact Illumina Technical Support.

TSO Comprehensive Assay Security Certificate

The TSO Comprehensive analysis module uses HTTPS to encrypt data connections to make sure that run data is private and secure. HTTPS is required for remote access of the instrument using a web browser from another machine in the same network. The TSO Comprehensive analysis module requires the installation of a TSO Comprehensive (JP) security certificate in addition to the NextSeq 550Dx instrument TSO Comprehensive analysis module security certificate.

NOTE If the TSO Comprehensive analysis module Security Patch is installed on a NextSeq 550Dx instrument, then remote access from the customer-supplied PC via web browser using HTTPS to the NextSeq 550Dx Local Run Manager web portal is disabled.

To install the TSO Comprehensive (JP) security certificate, do as follows.

1. Open TSO Comprehensive analysis module on your instrument.
2. Use the Tools menu to navigate to the Modules & Manifests screen.
3. Select **TSO Comp (JP) module**.
4. Download the TSO Comprehensive HTTPS Certificate.
5. Extract the contents of the zip file.
6. Right-click the BAT file and select **Run as administrator**.
7. Follow the prompts to finish the installation, and then restart your browser.

Regenerate Security Certificate

If the instrument name changed recently, or the instrument was moved to a new domain, you must regenerate the security certificate to regain access to the NextSeq 550Dx instrument and the TSO Comprehensive analysis module. For instructions on how to regenerate the TSO Comprehensive analysis module security certificate, refer to the *NextSeq 550Dx Instrument Site Prep Guide for Japan* (document # 1000000058597).

To regenerate the TSO Comprehensive (JP) security certificate, do as follows.

1. On the instrument, log in to the Windows operating system.
2. Using Windows File Explorer, navigate to the directory where the KB service is installed (for example, `C:\Illumina\Local Run Manager\Modules\TSOCompJP\[VersionNumber]\KBApiService\bin\Scripts`).
3. Right-click the BAT file and select **Run as administrator**.
4. Follow the prompts to finish the installation.
5. To connect to TSO Comprehensive analysis module from another device, download and install the regenerated certificate on the remote device.

Revision History

| Document | Date | Description of Change |
|--------------------------------|------------------|---|
| Document # 200049183 v02 | June 2025 | <ul style="list-style-type: none">• Corrected file types and file names.• Added a row to describe the Claims Package Version in the Analysis Details table• Removed the Tumor Type row from the table that describes what is included in the Low Depth Report header |
| Document # 200049183 v01 | March 2025 | <ul style="list-style-type: none">• Added the following:<ul style="list-style-type: none">• Clarification about how the KB is provided.• Equation to calculate copy number.• Instructions for regenerating the security certificate.• Clarifications about RNA Library QC fail results.• Disclaimers specifying TMB and MSI are out of scope. |
| Document # 200049183 v00 | February 2024 | Initial release. |

Technical Assistance

For technical assistance, contact Illumina Technical Support.

Website: www.illumina.com

Email: techsupport@illumina.com

Safety data sheets (SDSs)—Available on the Illumina website at support.illumina.com/sds.html.

Product documentation—Available for download from support.illumina.com.



Illumina, Inc.
5200 Illumina Way
San Diego, California 92122 U.S.A.
+1.800.809.ILMN (4566)
+1.858.202.4566 (outside North America)
techsupport@illumina.com
www.illumina.com

CE

IVD

EC REP



Illumina Netherlands B.V.
Steenoven 19
5626 DK Eindhoven
The Netherlands

FOR IN VITRO DIAGNOSTIC USE. FOR EXPORT ONLY.

© 2025 Illumina, Inc. All rights reserved.

illumina®