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Tabelle 2: Hybrid-Capture-Sondenpanels für Illumina Complete Long Read Prep with Enrichment, Human

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[Illumina Complete Long Read Prep with Enrichment, Human](#)

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Quellen

1. Mehio R, Ruehle M, Catreux S, et al. DRAGEN Wins at Precision- FDA Truth Challenge V2 Showcase Accuracy Gains from Alt-aware Mapping and Graph Reference Genomes. [illumina.com/science/genomics-research/dragen-wins-precisionfda-challenge-showcase-accuracy-gains.html](https://www.illumina.com/science/genomics-research/dragen-wins-precisionfda-challenge-showcase-accuracy-gains.html). Veröffentlicht 2020. Aufgerufen am 12. Januar 2023.
2. Illumina. Accuracy improvements in germline small variant calling with the DRAGEN Bio-IT Platform. [illumina.com/content/dam/illumina/gcs/assembled-assets/marketing-literature/dragen-v4-accuracy-app-note-m-gl-01016/dragen-v4-accuracy-app-note-m-gl-01016.pdf](https://www.illumina.com/content/dam/illumina/gcs/assembled-assets/marketing-literature/dragen-v4-accuracy-app-note-m-gl-01016/dragen-v4-accuracy-app-note-m-gl-01016.pdf). Veröffentlicht 2022. Aufgerufen am 1. Februar 2023.
3. Pacific Biosciences. Preparing DNA for PacBio HiFi sequencing—Extraction and quality control. [pacb.com/literature/technical-note-preparing-dna-for-pacbio-hifi-sequencing-extraction-and-quality-control/](https://www.pacb.com/literature/technical-note-preparing-dna-for-pacbio-hifi-sequencing-extraction-and-quality-control/). Veröffentlicht 2022. Aufgerufen am 12. Januar 2023.
4. Pacific Biosciences. Preparing whole genome and metagenome libraries using SMRTbell prep kit 3.0. [pacb.com/wp-content/uploads/Procedure-checklist-Preparing-whole-genome-and-metagenome-libraries-using-SMRTbell-prep-kit-3.0.pdf](https://www.pacb.com/wp-content/uploads/Procedure-checklist-Preparing-whole-genome-and-metagenome-libraries-using-SMRTbell-prep-kit-3.0.pdf). Veröffentlicht 2022. Aufgerufen am 5. Oktober 2023.
5. Oxford Nanopore Technologies. Ligation Sequencing Kit. <https://store.nanoporetech.com/us/productDetail/?id=ligation-sequencing-kit-v14>. Aufgerufen am 5. Oktober 2023.
6. Pacific Biosciences. Low Yield Troubleshooting Guide. [pacb.com/wp-content/uploads/Guide-Low-Yield-Troubleshooting.pdf](https://www.pacb.com/wp-content/uploads/Guide-Low-Yield-Troubleshooting.pdf). Veröffentlicht 2018. Aufgerufen am 12. Januar 2023.
7. Illumina. Illumina Complete Long Read Prep, Human – Datenblatt. [illumina.com/content/dam/illumina/gcs/assembled-assets/marketing-literature/illumina-long-read-prep-human-data-sheet-m-gl-01420/illumina-long-read-prep-data-sheet-m-gl-01420.pdf](https://www.illumina.com/content/dam/illumina/gcs/assembled-assets/marketing-literature/illumina-long-read-prep-human-data-sheet-m-gl-01420/illumina-long-read-prep-data-sheet-m-gl-01420.pdf). Veröffentlicht 2022. Aufgerufen am 22. September 2023.
8. Illumina. Comprehensive whole-genome sequencing with Illumina Complete Long Read Prep, Human – Technischer Hinweis. [illumina.com/content/dam/illumina/gcs/assembled-assets/marketing-literature/illumina-long-read-prep-human-tech-note-m-gl-01421/ilmn-long-read-hu-tech-note-m-gl-01421.pdf](https://www.illumina.com/content/dam/illumina/gcs/assembled-assets/marketing-literature/illumina-long-read-prep-human-tech-note-m-gl-01421/ilmn-long-read-hu-tech-note-m-gl-01421.pdf). Veröffentlicht 2022. Aufgerufen am 22. September 2023.
9. Bekritsky MA, Colombo C, Eberle MA. Identifying genomic regions with high quality single nucleotide variant calling. [illumina.com/science/genomics-research/articles/identifying-genomic-regions-with-high-quality-single-nucleotide-.html](https://www.illumina.com/science/genomics-research/articles/identifying-genomic-regions-with-high-quality-single-nucleotide-.html). Veröffentlicht 2021. Aufgerufen am 30. August 2023.
10. Wagner J, Olson ND, Harris L, et al. Curated variation benchmarks for challenging medically relevant autosomal genes. *Nat Biotechnol.* 2022;40(5):672-680. doi:10.1038/s41587-021-01158-1
11. PharmGKB. VIPs: Very Important Pharmacogenes. [pharmgkb.org/vips](https://www.pharmgkb.org/vips). Aufgerufen am 22. September 2023.
12. National Library of Medicine. GTR: Genetic Testing Registry. Precision HealthPGx Panel (25 Genes). ncbi.nlm.nih.gov/gtr/tests/593428/. Aktualisiert am 29. November 2022. Aufgerufen am 22. September 2023.
13. Pratt VM, Everts RE, Aggarwal P, et al. Characterization of 137 Genomic DNA Reference Materials for 28 Pharmacogenetic Genes: A GeT-RM Collaborative Project. *J Mol Diagn.* 2016;18(1):109-123. doi:10.1016/j.jmoldx.2015.08.005
14. Miller DT, Lee K, Abul-Husn NS, et al. ACMG SF v3.1 list for reporting of secondary findings in clinical exome and genome sequencing: A policy statement of the American College of Medical Genetics and Genomics (ACMG). *Genet Med.* 2022;24(7):1407-1414. doi:10.1016/j.gim.2022.04.006
15. Kulski JK, Suzuki S, Shiina T. Human leukocyte antigen super-locus: nexus of genomic supergenes, SNPs, indels, transcripts, and haplotypes. *Hum Genome Var.* 2022;9(1):49. doi:10.1038/s41439-022-00226-5
16. Roessler K. Illumina Complete Long Reads software analysis workflow for human WGS. [illumina.com/science/genomics-research/articles/complete-long-read-software-analysis.html](https://www.illumina.com/science/genomics-research/articles/complete-long-read-software-analysis.html). Veröffentlicht 2023. Aufgerufen am 22. September 2023.

Bestellinformationen		Bestellinformationen	
Produkt	Katalog-Nr.	Produkt	Katalog-Nr.
Illumina Complete Long Read Prep with Enrichment, Human (24 samples)	20113832	Illumina Unique Dual Indexes, LT (48 indexes, 48 samples)	20098166
Illumina Complete Long Read Prep with Enrichment, Human (96 samples)	20113833	Illumina DNA/RNA UD Indexes Set A, Tagmentation (96 indexes, 96 samples)	20091654
Illumina Complete Long Read Prep with Enrichment, Human Comprehensive Panel (24 samples)	20113834	Illumina DNA/RNA UD Indexes Set B, Tagmentation (96 indexes, 96 samples)	20091656
Illumina Complete Long Read Prep with Enrichment, Human Comprehensive Panel (96 samples)	20113835	Illumina DNA/RNA UD Indexes Set C, Tagmentation (96 indexes, 96 samples)	20091658
Illumina Human Comprehensive Panel (24 samples)	20113836	Illumina DNA/RNA UD Indexes Set D, Tagmentation (96 indexes, 96 samples)	20091660
Illumina Human Comprehensive Panel (96 samples)	20113837	Illumina Analytics – 1 iCredit	20042038
Illumina Custom Enrichment Panel v2 (32 µl, 120 bp)	20073953		
Illumina Custom Enrichment Panel v2 (384 µl, 120 bp)	20073952		
Illumina Custom Enrichment Panel v2 (1536 µl, 120 bp)	20111339		

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© 2024 Illumina, Inc. Alle Rechte vorbehalten. Alle Marken sind Eigentum von Illumina, Inc. bzw. der jeweiligen Inhaber. Spezifische Informationen zu Marken finden Sie unter www.illumina.com/company/legal.html.